alethia®

C. difficile

DNA Amplification Assay

DNA Amplification Assay for the Detection of
Cytotoxigenic *C. difficile* in Stool Specimens

REF 480050

IVD

Rx Only

INTENDED USE

The Alethia *C. difficile* DNA amplification assay, performed on the Alethia reader, is a qualitative in vitro diagnostic test for the direct detection of toxigenic *C. difficile* in human stool specimens from pediatric and adult patients suspected of having *Clostridium difficile*-associated disease (CDAD).

The Alethia *C. difficile* assay utilizes loop-mediated isothermal DNA amplification (LAMP)^{1,2} technology to detect the pathogenicity locus (PaLoc)³ of toxigenic *Clostridium difficile*. The *Clostridium difficile* PaLoc is a gene segment present in all known toxigenic *C. difficile* strains. The *C. difficile* PaLoc codes for both the Toxin A gene (*tcdA*) and the Toxin B gene (*tcdB*), has conserved border regions, and is found at the same site on the *C. difficile* genome for all toxigenic strains.³ The Alethia *C. difficile* assay detects the PaLoc by targeting a partial DNA fragment on the Toxin A gene. The *tcdA* target region was selected as an intact region remaining in all known A+B+ and A-B+ toxinotypes.

Alethia *C. difficile* is intended for use in hospital, reference or state laboratory settings. The device is not intended for nonlaboratory point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The Alethia *C. difficile* DNA molecular assay is based on loop-mediated amplification technology, which uses specifically designed primers to the PaLoc pathogenicity locus to provide for specific and continuous isothermal DNA amplification. A by-product of this amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. This presence of turbidity signifies a positive reaction while the absence of turbidity represents a negative reaction. The Alethia *C. difficile* assay contains primers that specifically amplify a 204 bp region of the conserved 5' sequence of the *tcdA* gene within the PaLoc of toxigenic *C. difficile* in diarrheal stool samples from patients suspected of having *C. difficile*-associated disease (CDAD). The results of the assay are determined using the Meridian Alethia Incubator / Reader.

The Alethia *C. difficile* DNA molecular assay kit includes a Sample Preparation Apparatus and a Test Device. The Sample Preparation Apparatus, used for specimen dilution and preparation, is filled with Phosphate Buffered Saline with formalin-treated *Staphylococcus aureus*. The Alethia *C. difficile* Test Device contains one dry reagent lysosphere in each of two chambers; a TEST Chamber with *C. difficile* specific primers and a CONTROL Chamber with *S. aureus* specific primers. Together, the formalin-treated *Staphylococcus aureus* in the Sample Preparation Apparatus and the *S. aureus* specific primers in the CONTROL Chamber lysosphere function as the Internal Control for the Alethia *C. difficile* assay. Each patient specimen is added to the Sample Preparation Apparatus via a Sample Collection Brush. The patient specimen is diluted in Phosphate Buffered Saline with formalin-treated *Staphylococcus aureus*, combining any *C. difficile* present in the patient specimen with *S. aureus*. The patient specimen and the *S. aureus* Control Target are processed in parallel through the remaining Specimen Preparation steps and through the Test Procedure. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The Control *S. aureus* DNA target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

BIOLOGICAL PRINCIPLES

Toxigenic *Clostridium difficile* is a major cause of antibiotic associated diarrhea and colitis and is the causative agent for virtually all cases of pseudomembranous colitis. Although about 2% of normal healthy adults are colonized with *C. difficile*, many patients acquire this organism through nosocomial infection. Exposure to most antibiotics is thought to allow proliferation of toxigenic *C. difficile* by disrupting the normal intestinal flora. Two large toxin proteins (TcdA [or toxin A] and TcdB [toxin B]) are thought to be the primary virulence factors of *C. difficile*. These toxins are encoded by two separate genes, named *tcdA* and *tcdB*, respectively. Together, with three additional genes, they form a 19.6 kb pathogenicity locus called PaLoc.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia *C. difficile* Sample Preparation Apparatus:** Sampling unit consisting of sample preparation chamber, dropper tip, cap and Sample Dilution Buffer (Phosphate Buffered Saline and formalin treated *Staphylococcus aureus*, with sodium azide (0.09%) as a preservative)
2. **Alethia Reaction Buffer:** Tris-buffered solution containing sodium azide (0.09%) as a preservative
3. **Alethia *C. difficile* Test Device:** Two separate chambers containing dry reagent lysospheres comprised of DNA polymerase, Deoxyribonucleotide Triphosphate (dNTPs), and either *C. difficile* specific primers (TEST Chamber) or *S. aureus* primers (CONTROL Chamber)
4. **Sample Collection Brushes**
5. **Alethia Heat Treatment Tubes**

MATERIALS PROVIDED SEPARATELY

Alethia External Control Kit, Catalog Number: 479920

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital Thermometer with Max/Min Temperature Memory (e.g. Traceable® Lollipop™ Waterproof/Shockproof Thermometer)
3. Vortex mixer
4. Interval timer
5. Micropipette capable of dispensing 50 µL
6. Micropipette capable of dispensing 200 µL
7. Alethia reader

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Follow Biosafety Level 2 and Good Laboratory practices during testing.⁴ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
3. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
4. Quality Control Programs for Molecular Testing Laboratories should be employed.⁵
5. Alethia Sample Dilution Buffer contains formalin inactivated-organisms and should be treated as potentially infectious.
6. The Alethia *C. difficile* Test Device contains lyophilized reagents. The protective pouch should not be open until ready to perform the assay.
7. The Alethia *C. difficile* Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-27 C.

SPECIMEN COLLECTION AND PREPARATION

Sample type: Unformed samples indicative of CDAD.

Human stool samples, unpreserved: Samples should be transported and stored at 2-8 C prior to testing. Samples should be tested as soon as possible, but may be held up to 24 hours at 21-27 C or 5 days at 2-8 C.

Human stool samples, preserved in Cary-Blair-based media: Samples should be transported and stored at 2-8 C prior to testing. Samples should be tested as soon as possible, but may be held up to 5 days at 2-8 C.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (21-27 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN PREPARATION

NOTE: Ensure that the Alethia instrument is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia reader Operator's Manual for further information regarding instrument set-up and operation.

1. Mix stool sample thoroughly.
2. Collect mixed sample specimen using Sample Collection Brush.
 - a. Liquid Stool: Immerse Sample Collection Brush completely into specimen.
 - b. Semi-solid Stool: Rotate Sample Collection Brush over specimen surface, lightly coating the surfaces of the brush bristles. Approximately one-half of the brush should be coated; over-collection of stool may lead to clogging of the sample collection device.
3. Add the Sample Collection Brush to the Alethia C. *difficile* Sample Preparation Apparatus containing Sample Diluent. Break-off the Sample Collection Brush at the 'Break' mark on the handle. Replace and secure the Sample Preparation Apparatus cap. For stool in Cary-Blair-based media, transfer 200 µL of specimen to the Alethia C. *difficile* Sample Preparation Apparatus and secure the cap. Vortex the Sample Preparation Apparatus for a minimum of 10 seconds. Specimen in Alethia C. *difficile* Sample Preparation Apparatus prior to heat-treatment may be held at 2-27 C for up to 24 hours prior to testing.
4. Remove the tip cap from the Sample Preparation Apparatus and squeeze five to ten drops of sample into a clean Alethia Heat Treatment Tube.
5. Repeat Sample Preparation Steps for all samples to be processed.
6. Heat the Alethia Heat Treatment Tubes containing mixed sample in a dry-bath/heat block at 95 +/- 5 C for 10 +/- 2 minutes. Use the digital thermometer and the interval timer to monitor temperature and time.
7. Remove the tubes from the dry-bath/heat block and vortex for approximately 10 seconds. Heat-treated samples may be held at 2-8 C for 4 hours OR may be frozen and stored at ≤ -20 C for 1 day prior to addition to Reaction Buffer. Heat-treated samples may be frozen and thawed once.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single Alethia reader run.

1. Transfer 50 µL of heat-treated sample to an appropriately labeled Alethia Reaction Buffer tube.
2. Vortex the Reaction Buffer Tube containing heat-treated sample for approximately 10 seconds.
3. Repeat steps 1 and 2 for all the samples to be analyzed before proceeding.
4. Remove 1 Alethia C. *difficile* Test Device from its protective pouch per sample. Carefully open the device by lifting the latch, holding the tubes such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
5. Using a micropipette, first transfer 50 µL from the Reaction Buffer tube containing the heat treated sample to the TEST chamber (left/White Bead), and then transfer 50 µL from the Reaction Buffer tube containing the heat treated sample to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
6. Gently mix the device by tapping or flicking, being careful to not introduce air bubbles or get the reaction mixture into the device lid. Carefully examine the reaction tubes for air bubbles or liquid in the lid. If air bubbles or liquid is noted, tap the device on the bench top and repeat the visual inspection.
7. Insert the Alethia Test Device into the Alethia reader and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains toxigenic <i>C. difficile</i> strain with the pathogen locus (PaLoc).
	NEGATIVE	No toxigenic <i>C. difficile</i> detected.
	INVALID	No reportable result. Repeat the test using the original stool sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original stool samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original stool samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Each device contains an internal control well that controls for amplification inhibition, assay reagents and sample processing effectiveness.
2. The heat treatment step is monitored with an external thermometer and interval timer. Use the Max/Min temperature memory of the thermometer to ensure that a temperature of 95 +/- 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 +/- 2 minutes.
3. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
4. Alethia C. *difficile* External Control Reagents are supplied separately (Catalog 479920). It is recommended that the reactivity of each new lot and each new shipment of Alethia C. *difficile* be verified on receipt or before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia C. *difficile* test kit should not be used in patient testing if the external controls do not produce the correct results.
5. A separate device must be used for each external control reagent.
6. DO NOT heat-treat the Positive or Negative Control samples.

EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors, including patient population, type of institution and epidemiology. The incidence of CDAD in patients suspected of having antibiotic-associated disease is 15-20% although different facilities may find positive rates above or below this range.⁶ The incidence of toxinogenic *C. difficile* during the 2010 period of this study was 15.5% and 16.1%.

LIMITATIONS OF THE PROCEDURE

1. Assay does not distinguish between viable and non-viable organisms.
2. Test results are to be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
3. This test is for use with unfrozen stool samples, preserved or unpreserved, only. Performance characteristic of other clinical specimens types have not been established.
4. Two distinct groups have been identified that can harbor *C. difficile* asymptomatically at very high rates. Colonization at rates up to 50% and higher have been reported in infants and rates up to 32% in cystic fibrosis patients.^{7, 8}
5. This assay does not identify antimicrobial susceptibility.
6. This test detects but does not differentiate the NAP1 (Ribotype 027) strain from other toxigenic strains of *C. difficile*.
7. The detection of bacterial nucleic acid is dependent on proper specimen collection, handling (including transportation and storage) and preparation (dilution and heat-treatment). Failure to follow instructions for collection, handling and preparation may cause incorrect results.

SPECIFIC PERFORMANCE CHARACTERISTICS

Alethia *C. difficile* was evaluated in 2010 by four independent clinical test sites located in the Midwestern and Southern regions of the United States and the manufacturer. Two individual studies were performed to establish performance characteristics for the assay. Clinical study information is described below.

(1) Patients 2 Years of Age and Above

A total of 697 qualified samples from patients 2 years of age and above were evaluated. Samples were collected from 274 (39.3%) males and 419 (60.1%) females. In the case of 4 (0.6%) of the patients, gender was not known. The age groups of patients range from 2 years of age to 96 years of age. No differences in test performance were observed based on patient age or gender. The overall positivity rate for the study was 15.5%. Table 1 shows overall performance data for the study. Table 2 shows assay performance by clinical site.

Table 1. Performance data (Patients 2 years of age and above)

Cytotoxic bacterial culture	Alethia <i>C. difficile</i>			
	Positive	Negative	Invalid***	Total
Positive	99	5**	4	108
Negative	27*	546	16	589
Total	126	551	20	697
			95% CI	
Sensitivity	99/104	95.2%	89.2 - 97.9%	
Specificity	546/573	95.3%	93.2 - 96.7%	
Correlation	645/677	95.3%	93.4 - 96.6%	
Invalid Rate	20/697	2.9%	N/A	

* 15/27 false positive results were positive by another FDA cleared molecular assay. Of the remaining 12 false positive results, 8 were positive by a FDA cleared assay for the detection of GDH.

** 2/5 false negative results were negative by another FDA cleared molecular assay.

*** Invalid results were obtained for 20/697 (2.9%) samples tested. Eleven (1.6%) of the invalids observed were categorized as Assay Invalids, indicative of improper sample preparation, reagent failure, instrument failure or internal control failure. One of the 11 specimens remained invalid after repeat testing from the original sample.

Table 2. Performance characteristics by site (Patients 2 years of age and above)

Site	Positive Samples			Negative Samples		
	Alethia / Cytotoxic bacterial culture	Sensitivity %	95% CI	Alethia / Cytotoxic bacterial culture	Specificity %	95% CI
Total	99/104**	95.2%	89.2 – 97.9%	546/573*	95.3%	93.2 – 96.7%
Site 1	4/5	80.0%	37.6 – 96.4%	58/60	97.6%	88.6 – 99.1%
Site 2	12/12	100%	75.7 – 100%	62/67	92.5%	83.7 – 96.8%
Site 3	20/20	100%	83.9 – 100%	87/92	94.6%	87.9 – 97.7%
Site 4	8/8	100%	67.6 – 100%	36/39	92.3%	79.7 – 97.3%
Site 5	55/59	93.2%	83.8 – 97.3%	303/315	96.2%	93.5 – 97.8%

(2) Patients Less than 2 Years of Age

A total of 193 qualified samples from patients 0 months to 24 months were evaluated. Samples were collected from 103 (53.4%) males and 90 (46.6%) females. No differences in test performance were observed based on age or gender. The overall positivity rate for the study was 16.2%. Table 3 shows overall performance data for the study. Table 4 shows assay performance by clinical site.

Table 3. Performance data (Patients less than 2 years of age)

Cytotoxic bacterial culture	Alethia <i>C. difficile</i>			
	Positive	Negative	Invalid***	Total
Positive	28	2**	1	31
Negative	6*	156	0	162
Total	34	158	1	193
			95% CI	
Sensitivity	28/30	93.3%	78.7 - 98.2%	
Specificity	156/162	96.3%	92.2 - 98.3%	
Correlation	184/192	95.8%	91.3 - 97.5%	
Invalid Rate	1/193	0.5%	N/A	

* 3/6 false positive results were positive by another FDA cleared molecular assay. Of the remaining 3 false positive results, all were positive by a FDA cleared assay for the detection of GDH.

** 1/2 false negative results were negative by another FDA cleared molecular assay.

*** Invalid results were obtained for 1/193 (0.5%) samples tested. The invalid observed was categorized as an Assay Invalid, indicative of improper sample preparation, reagent failure, instrument failure or internal control failure. The specimen remained invalid after repeat testing from the original sample.

Table 4. Performance characteristics by site (Patients less than 2 years of age)

Site	Positive Samples			Negative Samples		
	Alethia / Cytotoxic bacterial culture	Sensitivity %	95% CI	Alethia / Cytotoxic bacterial culture	Specificity %	95% CI
Total	28/30**	93.3%	78.7 - 98.2%	156/162	96.3%	92.2 - 98.3%
Site 1	8/8	100%	67.6 - 100%	48/49	98.0%	89.3 - 99.6%
Site 2	20/22	90.9%	72.2 - 97.5%	105/109	96.3%	90.9 - 98.6%
Site 4	0/0	N/A	N/A	2/3	66.7%	20.8 - 93.9%
Site 5	0/0	N/A	N/A	1/1	100%	20.7 - 100%

Table 5. Overall results by patient age

Patient age	Positive Samples			Negative Samples		
	Alethia / Toxigenic culture	Sensitivity %	95% CI	Alethia / Toxigenic culture	Specificity %	95% CI
< 2 years	28/30	93.3%	78.7 - 98.2%	156/162	96.3%	92.2 - 98.3%
≥ 2 - 12 years	10/11	90.9%	62.3 - 98.4%	75/79	94.9%	87.7 - 98.0%
> 12 to 21 years	5/5	100%	56.6 - 100%	53/56	94.6%	85.4 - 98.2%
> 21 years	83/87	95.4%	88.8 - 98.2%	417/437	95.4%	93.0 - 97.0%
Age Unknown	1/1	100%	20.7 - 100%	1/1	100%	20.7 - 100%

ANALYTICAL SENSITIVITY

The analytical sensitivity of this assay for *C. difficile* was based on 20 replicates for each measurand and with a stated probability (eg, 95% where 19/20 replicates are positive) of obtaining positive responses at the following levels of the measurands:

Strain ID	Toxinotype	Phenotype	LoD/Test
VPI 10463	0	A+/B+	4 CFU/test
2007431	III (NAP1)	A+/B+	32 CFU/test
CFI	VIII	A-/B+	64 CFU/test
2006240	V (NAP7)	A+/B+	32 CFU/test
BI8	III	A+/B+	64 CFU/test
2007858	IX/XXIII	A+/B+	32 CFU/test
8864	X	A-/B+	64 CFU/test

ASSAY REACTIVITY

The following *C. difficile* stock cultures from different sources were tested and produced positive reactions at 64 CFU/test with Alethia *C. difficile*:

Type 0: Strains 10463, 2004111, 2004205, 2005070, 2005257, 2008029, 2008162, 2008341, 2008351, 2009066, 2009099, B1, G1, J7, K12, Y1; Type III: 2004052, 2004118, 2007431, BI17, BI8; Type V: 2005325, 2006240, 2008188, 2009018, 2009065, BK6; Type VIII: 43598, 2008016, CFI; Type X 8864; Type XII 2007435; Type IX/XXIII 2007858; Unknown Type 2009132, 2009155, 2009277.

REPRODUCIBILITY

Blind coded panels of 10 samples were supplied to three independent laboratories for precision studies. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured at the assay limit of detection (n = 3) and just below the limit of blank (ie, high negative sample, n = 3). The panels also included uncharacterized positive (n = 2) and negative (n = 2) samples. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of Alethia *C. difficile* were used in this study. The results are given in the table below:

Sample Type	Site 1		Site 2		Site 3		Total	
	Percent agreement		Percent agreement		Percent agreement		Percent agreement	
Negative	20/20	100%	20/20	100%	19/19****	100%	59/59	100%
High Negative	25/30	83%	29/30	97%	28/30	93%	82/90	91%
Low Positive	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Positive	20/20	100%	20/20	100%	20/20	100%	60/60	100%

**** 1 specimen generated an instrument invalid test result.

CROSSREACTIVITY STUDIES

Crossreactivity studies were performed with positive and negative stool specimens inoculated with bacterial or fungal organisms to a final concentration of 1.2×10^8 or virus at a minimum of $1 \times 10^{5.06}$ TCID₅₀/mL. None of the following organisms in stool reacted with Alethia *C. difficile*:

Aeromonas hydrophila, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter freundii*, *Clostridium sordellii*, *Clostridium perfringens*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* O157:H7, *Escherichia fergusonii*, *Escherichia hermannii*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Listeria monocytogenes*, *Peptostreptococcus anaerobius*, *Plesiomonas shigelloides*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Salmonella* Groups B-E, *Serratia liquefaciens*, *Serratia marcescens*, *Shigella boydii*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*, Adenovirus Types 40 and 41, Coxsackievirus, Echovirus, Rotavirus.

TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluents concentrations, do not interfere with test results in the final concentrations listed: Barium sulfate (5 mg/mL), fecal fat (equivalent to 2.65 mg stearic plus 1.3 mg palmitic acids per mL), hemoglobin (as methemoglobin) (3.2 mg/mL), IgA (5 mg/mL), Imodium AD® (0.00667 mg/mL), Kaopectate® (0.87 mg/mL), Metronidazole (12.5 mg/mL), mucin (3.33 mg/mL) Mylanta® (4.2 mg/mL), Pepto-Bismol® (0.87 mg/mL), Prilosec® (0.5 mg/mL), Tagamet® (0.5 mg/mL), TUMS® (0.5 mg/mL), Vancomycin (12.5 mg/mL), white blood cells (5% V/V), whole blood (5% V/V).

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REFERENCES

1. Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes 2002;16:223-29.
2. Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. J Biochem Biophys 2004;59:145-47.
3. Rupnik M. Heterogeneity of large clostridial toxins: importance of Clostridium difficile toxinotypes. Microbiol Rev 2008;32:541-55.
4. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
5. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.
6. Wilkins TD, Lysterly DM. Clostridium difficile testing: after 20 years, still challenging. J Clin Microbiol 2003;41:531-4.
7. Johnson S, Gerding DN. Clostridium-difficile-associated diarrhea. Clin Infect Dis 1998;26:1027-34.
8. Lysterly DM, Krivan HC, Wilkins TD. Clostridium difficile: its disease and toxins. Clin Microbiol Rev 1988;1:1-18.

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
EC REP	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: - For Italian Customers: ordini@meridianbioscience.com - For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
AUSTRIAN SPONSOR	Emergo Australia Level 20, Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia
CH REP	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	Meridian products carrying the European Conformity (CE) mark fulfill the requirements of Directive 98/79/EC or the Regulation 2017/746 on in-vitro diagnostic medical devices / I prodotti Meridian recanti il marchio di Conformità Europea (CE) soddisfano i requisiti della Direttiva 98/79/CE o del Regolamento 2017/746 sui dispositivi medico-diagnostici in vitro / Les produits Meridian portent la marque de Conformité européenne (CE) sont conformes aux exigences de la Directive 98/79/CE ou du Règlement 2017/746 relatifs aux dispositifs médicaux de diagnostic in vitro / Los productos de Meridian que llevan la marca de conformidad europea (CE) cumplen los requisitos de la Directiva 98/79/CE o del Reglamento 2017/746 sobre los productos sanitarios para diagnóstico in vitro / Produkte von Meridian mit der CE-Kennzeichnung erfüllen die Anforderungen der EU-Richtlinie 98/79/EG bzw. der Verordnung (EU) 2017/746 über In-vitro-Diagnostika	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
		STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Stérilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Cautions, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Stérilisation par oxyde d'éthylène / Esterilizado por óxido de etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contenu suffisant pour <n> tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femminile / De sexe féminin / Hembra / Frau		Recycle – do not dispose of as general waste / Riciclare – non eliminare come rifiuto generico / Recycler – ne pas jeter dans une poubelle / Recicle – no deseché como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung	TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	HT TUBE	Heat Treatment Tube / Provetta per il Trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrchen zur Hitzebearbeitung
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Eviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden		HOT SURFACE: Keep hands away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	IPX-0	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	MEDIA	Media / Terreno di trasporto / Milieux / Medio / Medium
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrchen mit Schnappverschluss	COL	Sample Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampón de l'échantillon / Tampón de muestra / Probenpuffer	PRE REAG	Pretreatment Reagent / Reagente di Pretrattamento / Réactif de prétraitement / Reactivo de pretratamiento / Reagenz für die Vorbehandlung
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung
TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß	CH REP	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.