FOR EXPORT ONLY. Not for sale in the USA.

alethia[®]

Chlamydia

DNA Amplification Assay

DNA Amplification Assay for the Detection of *Chlamydia trachomatis*

REF 480450

IVD

INTENDED USE

The Alethia Chlamydia (CT) assay, performed on the Alethia Incubator/Reader, is an in vitro diagnostic device for the detection of *Chlamydia trachomatis* in physician-collected (or other trained medical professional) endocervical swabs, patient-collected vaginal swab specimens from clinical settings, male urine specimens and female urine specimens.

The Alethia CT assay utilizes loop-mediated isothermal DNA amplification (LAMP) technology to detect *Chlamydia trachomatis* by targeting a segment of the *Chlamydia trachomatis* genome. Results from the Alethia CT assay can be used as an aid in the diagnosis of chlamydial urogenital infection in asymptomatic and symptomatic individuals. This assay is not intended to monitor treatment for chlamydial infections.

The Alethia CT assay is intended for use in hospital, reference or state laboratory settings. The device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The Alethia Chlamydia molecular assay is based on loop-mediated amplification (LAMP) technology.¹ The assay targets a 240 base pair (bp) sequence of *Chlamydia trachomatis* cryptic plasmid. The target DNA sequence is found in a region of the plasmid which codes for an unnamed protein product and is found in all 33 strains of *C. trachomatis* in which the cryptic plasmid sequence is reported.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia Chlamydia kit includes Alethia Chlamydia/Gonorrhea Assay Control, Alethia Chlamydia Test Devices and Mineral Oil. The Alethia Assay Control, used for specimen dilution and preparation, is a buffered solution containing formalin-treated *Escherichia coli* cells harboring *Staphylococcus aureus* DNA. The Alethia Chlamydia Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with Chlamydia-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA from the Control Reagent and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, each patient specimen is diluted with the Control Reagent containing the control *S. aureus* target. Addition of *S. aureus* DNA to the patient specimen allows for parallel processing of target DNA and Control DNA through sample preparation, amplification and detection. The Internal Control monitors DNA extraction, amplification inhibition, assay reagent performance and sample processing effectiveness. The control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The Alethia Incubator/Reader monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{initial}, S_i) and at the assay Run End (Signal_{final}, S_f). The Alethia reader calculates the change in light transmission between Run End and Run Start (S_f:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_f:S_i ratios less than 82% are reported as 'POSITIVE'; TEST chamber S_f:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_f:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_f:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Chlamydia trachomatis is one of the most common bacterial sexually transmitted infections worldwide. The World Health Organization estimates that 92 million cases of *C. trachomatis* infection occur annually.² Chlamydia remains the most commonly reported infectious sexually transmitted disease in the United States, with more than 1.3 million cases reported in 2010.³

Chlamydiae are nonmotile, gram-negative, obligate intracellular bacteria. The *C. trachomatis* species is comprised of fifteen serovars (A, B, Ba, C, D, E, F, G, H, I, J, K, L1, L2 and L3), categorized into 3 groups, each associated with a different disease state. Serovars A-C are responsible for ocular disease; Serovars D through K are the major cause of genital chlamydial infections in men and women. Serovars L1, L2, and L3 are responsible for lymphogranuloma venereum. *C. trachomatis* infection causes cervicitis, urethritis, salpingitis, proctitis, and endometritis in women, and urethritis, epididymitis and proctitis in men. Asymptomatic infection is common in both males and females. Left untreated, chlamydia can cause serious health problems in women such as Pelvic Inflammatory Disease (PID), ectopic pregnancy, and infertility.⁴ The organism can be transmitted in the birth canal; children born to infected mothers are at significantly higher risk for neonatal ophthalmia and pneumonia.⁴ The United States Centers for Disease Control and Prevention recommend annual chlamydia screening for all sexually active women ≤ 25 years due to the risks associated with chlamydial infection.⁴

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia Chlamydia/Gonorrhea Assay Control (Assay Control):** Phosphate-buffered solution containing a blue dye and formalin-treated *E. coli* cells harboring plasmid containing a segment of the *S. aureus* genome with sodium azide (0.09%) as a preservative. Packaged in bottle with dropper tip.
2. **Alethia Chlamydia Test Device:** Two-chambered device packaged in a foil pouch containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either Chlamydia-specific primers (TEST Chamber) or control primers (CONTROL Chamber).
3. **Mineral Oil**
4. **Transfer Pipettes**

MATERIALS PROVIDED SEPARATELY

M-prep™ DNA Sample Preparation Kit, Meridian Bioscience, Inc. Catalog Number: 281110

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Sterile, plastic, preservative-free, urine specimen collection cups
4. MSwab™ Collection, Preservation and Transport System for Molecular and Culture Applications, Copan Flock Technologies Srl, Catalog Number 6E028N.MER
5. NATtrol™ *Chlamydia trachomatis* External Run Controls (Strain D-UW-3), Zeptomatrix, Catalog Number NATCT(D-UW3)-ERCL
6. Molecular Biology Grade Tris-EDTA (TE) Buffer
7. Nuclease-free Microcentrifuge Tubes, 2 mL

EQUIPMENT NOT PROVIDED

1. Microcentrifuge capable of 2000 ± 250 rcf (2000 ± 250 x g)
2. Alethia Incubator/Reader, Meridian Bioscience, Inc. Catalog Number: 610189
3. Micropipette(s) capable of dispensing 50 µL and 250 µL
4. Interval Timer

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange Assay Control Reagent or Test Devices between lots. Mineral Oil is interchangeable provided it is within assigned expiration dates when used.
3. Follow Biosafety Level 2 and Good Laboratory practices during testing.⁵ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
5. Quality Control Programs for Molecular Testing Laboratories should be employed.⁵
6. The Alethia Chlamydia Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
7. The Alethia Chlamydia Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS


There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

PROCEDURAL NOTES

The transfer pipette included in the Alethia Chlamydia kit is individually wrapped and capable of dispensing up to 3 mL. The pipette is diagrammed below.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-30 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN COLLECTION AND PREPARATION

Sample types: Physician-collected endocervical swabs, patient-collected vaginal swabs (from clinical settings), male urine and female urine

Endocervical and Vaginal Swab Specimen Collection

Endocervical and vaginal swab specimens should be collected by trained medical professionals, in accordance with institutional guidelines.

Vaginal swab specimens may be self-collected by the patient when instruction and supervision is provided by a trained medical professional. Institutional guidelines should be followed.

Urine Specimen Collection

Urine specimens must be collected in a sterile, plastic, preservative-free, specimen collection cup. The following steps should be taken to ensure specimen quality:

1. Do not let patient urinate for at least 1 hour prior to specimen collection.
2. Have patient collect the first 10-30 mL of voided urine (the first part of the stream, not midstream) directly into the collection cup.
3. Cap and label the collection cup with patient and sample collection information.

Collected swab and urine specimens should be stored refrigerated (2-8 C) after collection. Refrigerated specimens may be held for up to 7 days prior to testing. Do not freeze specimens.

SPECIMEN PREPARATION

NOTE: Turn on the Alethia Incubator/Reader. Ensure that the required performance verifications have been completed before initiation of SPECIMEN PREPARATION. Refer to the Alethia reader Operator's Manual for further information regarding instrument set-up and operation.

Refer to the *M*-prep DNA Sample Preparation Kit package insert for information regarding column preparation, orientation and use. Use 1 *M*-prep column per sample. Remove the column top cap, followed by the bottom twist-off tip, and place the column tip into an *M*-prep tube to drain. Drained *M*-prep columns should be used within 1 hour.

Ensure specimens are at room temperature (19-30 C) before specimen preparation.

Swab Specimen Preparation

1. Swirl or invert the collected swab in medium to mix. Open the medium container carefully.

2. Using a transfer pipette, transfer 1.5 mL of the swab medium into a 2.0 mL screw-top Microcentrifuge *M*-prep tube. Discard the transfer pipette.
3. Briefly vortex the Alethia Chlamydia/Gonorrhea Assay Control. Add 1 drop of the Assay Control to the *M*-prep tube containing the sample medium. **The prepared sample will have a light blue appearance after addition of the Assay Control reagent.**
4. Cap the *M*-prep tube containing the sample.

Urine Specimen Preparation

1. Swirl the urine collection cup to mix. Open the urine container carefully.
2. Using a transfer pipette, transfer 1.5 mL of urine into a 2.0 mL screw-top microcentrifuge *M*-prep tube. Discard the transfer pipette.
3. Briefly vortex the Alethia Chlamydia/Gonorrhea Assay Control. Add 1 drop of the Assay Control to the *M*-prep tube containing the sample medium. **The prepared sample will have a light blue appearance after addition of the Assay Control reagent.**
4. Cap the *M*-prep tube containing the sample.

External Positive Control Preparation

1. Vortex NATrol™ *Chlamydia trachomatis* External Run Control.
2. Transfer 1.5 mL TE Buffer to a 2.0 mL nuclease-free microcentrifuge tube.
3. Using a micropipette, add 250 µL External Run Control to the tube containing TE Buffer.
4. Briefly vortex the Alethia Chlamydia/Gonorrhea Assay Control. Add 1 drop of the Assay Control to the *M*-prep tube containing the sample medium.
5. Cap the microcentrifuge tube containing External Control.

External Negative Control Preparation

1. Transfer 1.75 mL TE Buffer to a 2.0 mL nuclease-free microcentrifuge tube.
2. Briefly vortex the Alethia Chlamydia/Gonorrhea Assay Control. Add 1 drop of the Assay Control to the *M*-prep tube containing the sample medium.
3. Cap the microcentrifuge tube containing External Control.

Sample Processing

NOTE: The amount of prepared sample collected at the end of the Sample Processing steps should be enough to complete 2 Alethia tests.

1. Centrifuge each *M*-prep tube containing sample for 3 minutes at 2000 ± 250 rcf (2000 ± 250 x g).
2. Remove the *M*-prep tube from the microcentrifuge. Manually decant the supernatant fluid without disturbing any pellet formed during centrifugation.
3. Using a micropipette, add 250 µL of *M*-prep Buffer I to the *M*-prep tube. Discard the pipette tip. Vortex 10 seconds to mix.
4. Using a micropipette set to 250 µL, transfer the entire prepared sample to the top of an appropriately labeled and prepared *M*-prep Column. Wait approximately 2 minutes, or until the sample has been absorbed by column and flow stops. This step should take no longer than 30 minutes.
5. Using a micropipette, add 250 µL of *M*-prep Buffer II to the top of the *M*-prep column. Discard the pipette tip. **The column will have a red appearance after the addition of *M*-prep Buffer II.** Wait approximately 2 minutes, or until the red-colored buffer is absorbed by the column and flow stops. This step should take no longer than 30 minutes.
6. Remove the last drop of liquid from the column tip with the *M*-prep Tube. Discard the tube.
7. Place a clean *M*-prep Tube under the *M*-prep Column. Using a micropipette, add 250 µL of *M*-prep Buffer III to the top of the *M*-prep Column. Discard the pipette tip. Wait approximately 2 minutes or until flow stops. This step should take no longer than 30 minutes.
8. Remove the last drop of liquid from the column tip with the *M*-prep Tube. Label the tube with sample identification information and proceed to Test Procedure.

TEST PROCEDURE

NOTE: A maximum of 10 test devices can be processed in a single Alethia reader run.

1. Remove 1 Alethia Chlamydia Test Device from its protective pouch per sample. Carefully open the device by lifting the latch, holding the device upright so that the lyophilized reagents do not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
2. Transfer 50 µL of sample to the TEST chamber (Left/White Bead) and CONTROL chamber (Right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette.
3. Add 1 drop of Mineral Oil to both the TEST chamber and CONTROL chamber. Close the Alethia Test Device and fasten the latch securely.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the Test Device for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
5. Insert the Alethia Test Device into the Alethia reader and initiate run. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Chlamydia trachomatis</i> target DNA.
	NEGATIVE	No <i>Chlamydia trachomatis</i> DNA detected.
	INVALID	No reportable result. Repeat the test using the original patient specimen. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- 1. Each device contains an internal control that controls for amplification inhibition, assay reagents and sample processing effectiveness.
- 2. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- 3. It is recommended that the reactivity of each new lot and each new shipment of Alethia Chlamydia be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia Chlamydia test kit should not be used in patient testing if the external controls do not produce the correct results.
- 4. A separate Test Device must be used for each external control reagent.

EXPECTED VALUES

The incidence of *Chlamydia trachomatis* as detected by the Alethia Chlamydia Assay in prospectively collected, non-selected specimens during the period of clinical studies was 2.8% (3/107) with endocervical swabs, 3.9% (4/102) with vaginal swabs, and 16.1% (73/454) in urine.

LIMITATIONS OF THE PROCEDURE

- 1. This product can only be used with the Alethia instrument.
- 2. The Alethia Chlamydia assay is a qualitative assay and does not provide quantitative values or information about organism load.
- 3. Performance of the Alethia Chlamydia assay with endocervical and vaginal swabs has been established using the Copan MSwab™ Collection, Preservation and Transport System for Molecular and Culture Applications, Catalog 6E028N.MER. The use of alternate swab collection systems may lead to improper results and is not recommended.
- 4. The detection of nucleic acids is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
- 5. Organism nucleic acid may persist *in vivo*, independent of organism viability. The Alethia Chlamydia assay does not distinguish between viable and nonviable organisms.
- 6. As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Alethia Chlamydia DNA Amplification Assay was evaluated between July and September 2014. Clinical sample collection was carried out by three independent facilities representing geographically distinct regions within the United States. Physician-collected endocervical swabs, patient-collected vaginal swabs, female urine, and male urine specimens were collected prospectively from symptomatic and asymptomatic patients. The study was supplemented with left-over clinical urine specimens taken from patients for whom chlamydia and/or gonorrhea testing was ordered by a physician. Clinical specimens were tested with Alethia Chlamydia by the manufacturer.

The performance of Alethia Chlamydia was compared to an established, commercially available PCR assay for the detection of *Chlamydia trachomatis* and with a comparator method by an independent Reference Laboratory. Clinical study information is described below.

Table 1 Performance Characteristics - *C. trachomatis*

Comparator	Alethia Chlamydia All Samples Combined, All Patients Female Endocervical Swabs, Female Vaginal Swabs, Female Urine, Male Urine			
	Positive	Negative	Invalid	No Flow
Positive	67	2	4	0
Negative	12	550	12	6
Not Tested	1	9	0	0
Total	80	561	16	6
95% CI				
Sensitivity	67/69	97.1%	90.0 – 99.2%	
Specificity	550/562	97.9%	96.3 – 98.8%	
Correlation	617/631	97.8%	96.3 – 98.7%	
Invalid Rate	16/663	2.4%	1.5 – 3.9%	

Sixteen specimens produced invalid results, five swabs and 11 urines. Six specimens did not flow through the *M-Prep* column and did not produce Alethia results.

Samples used in the study were taken from patients for whom chlamydia and/or gonorrhea testing was ordered by a physician. One hundred seventeen (25.0%) patients were considered asymptomatic and 55 (11.8%) were considered symptomatic. The remaining 295 patients were not classified.

Age information was known for 466 of the patients from whom samples were tested. Patient age ranged from 14 years to 66 years. Patients less than 14 years or age were excluded from the study. Two hundred ninety nine (64.0%) patients were between 14 and up to 21 years of age; 166 (35.5%) were between 21 and up to 65 years old; and 1 (0.2%) was above 65 but below 90 years. Age for 1 (0.2%) patient was unknown. There were no noted performance differences based on age.

The study population included 277 (59.3%) female and 190 (40.7%) male patients. There is no expectation that the Alethia Chlamydia assay performance characteristics are influenced by patient gender.

ANALYTICAL SENSITIVITY

The analytical sensitivity or Limit of Detect for the Alethia Chlamydia assay was determined for one *C. trachomatis* strain, representing serovar D. Limit of Detection was determined using 60 replicates and a stated probability (eg. 95%, where 57/60 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below.

Serovar	Strain	IFU/mL
D	D-UW3r	2.5

ASSAY REACTIVITY

The following *C. trachomatis* strains were tested and produced a positive reaction with Alethia Chlamydia at the indicated concentrations:

Serovar	Strain ID	Concentration (IFU/mL)
A	Har 13	2.5
B	Har 36	0.625
Ba	Apache 2	0.625
C	Har 32	0.625
E	Bour	10.0
F	IC-Cal3	0.625
G	UW 57	0.625
H	580	0.625
I	UW12	0.625
J	UW36	1.25
K	UW31	0.625

REPRODUCIBILITY

Reproducibility studies were carried out by three laboratories. Blind-coded panels of 10 samples were supplied to participating laboratories. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured as moderate positive samples (15 IFU/mL), low positive samples (5 IFU/mL); and high negative samples (0.078 IFU/mL). The panel also included one negative sample, positive control and negative control. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of Alethia Chlamydia and nine Alethia readers were used in this study. Positive and Negative Controls were tested each day of testing. The results are provided in the table below:

Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive	29/30	96.7%	29/30	96.7%	28/30	93.3%	86/90	95.6%
Low Positive	30/30	100.0%	29/30	96.7%	26/30	86.7%	85/90	94.4%
High Negative	30/30	100.0%	30/30	100.0%	29/30	96.7%	89/90	98.9%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Negative Control	10/10	100.0%	10/11	90.9%	10/10	100.0%	30/31	96.8%
Positive Control	10/10	100.0%	10/11	90.9%	10/10	100.0%	30/31	96.8%

CROSSREACTIVITY

Crossreactivity studies were performed with positive and negative urine specimens inoculated with microorganisms to a final concentration of 1.2×10^6 CFU/mL. None of the following organisms reacted with Alethia Chlamydia: *Acinetobacter lwoffii*, *Alcaligenes faecalis*, *Bacteroides fragilis*, *Candida albicans*, *Candida tropicalis*, *Corynebacterium xerosis*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Listeria monocytogenes*, *Morganella morganii*, *Neisseria cinerea*, *Neisseria denitrificans*, *Neisseria elongata*, *Neisseria flava*, *Neisseria flavescens*, *Neisseria lactamica*, *Neisseria meningitidis*, *Neisseria mucosa*, *Neisseria perflava*, *Neisseria polysaccharea*, *Neisseria sicca*, and *Neisseria subflava*.

TESTS FOR INTERFERING SUBSTANCES

The following chemical substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results: Acetylsalicylic Acid (4 mg/mL), Acetaminophen (3.2 mg/mL), and Feminine Body Powder (2.5% w/v).

The following biological substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results: Acidic Urine [pH 4, (100% v/v)], Alkaline Urine [pH 9, (100% v/v)], Beta Estradiol (0.07 mg/mL), Human DNA (3.333 ng/mL), Leukocytes (1.04×10^6 cells/mL), Mucin [bovine submaxillary gland type I-S (0.2% w/v)], Progesterone (7 mg/mL), Seminal Fluid (5.0% v/v), Whole Blood (0.3% v/v).

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- REFERENCES**
1. Nagamine, K., Hase, T., Notoni, T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes 2002;16:223-29.
 2. World Health Organization. 2005. Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. World Health Organization, Geneva, Switzerland.
 3. Centers for Disease Control and Prevention. Sexually Transmitted Diseases Surveillance 2010 Atlanta: U.S Department of Health and Human Services; 2011.
 4. Centers for Disease Control and Prevention. Sexually Transmitted Diseases Treatment Guidelines, 2012. MMWR 2010; 59(No. RR_12): [44-55].
 5. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

0459

SN11031

REV. 08/21

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S.r.L.
Via dell'Industria, 7
20035 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bnl@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe B.V.
Zekeringstraat 17 A
1014BM Amsterdam
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: Info.bnl@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositif pour la préparation du campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdüngungspuffer befindet
		STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Sterilisation par irradiation aux rayons gamma / Esterilizado por radiación gamma / Sterilisation durch Gammastrahlen
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Caution, consult accompanying documents / Attention, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Stérilisation par oxyde d'éthylène / Esterilizado por óxido de etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contenu suffisant pour <n> tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction / tamponee / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femminile / De sexe féminin / Hembra / Frau		Recycle – do not dispose of as general waste / Riciclare – non eliminare come rifiuto generico / Recycle – ne pas éliminer dans une poubelle / Recycle – no deshecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung		Heat Treatment Tube / Provetta per il Trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrenchen zur Hitzebehandlung
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät		HOT SURFACE: Keep hands Away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Éviter l'exposition au faisceau / Radiación Laser: Evite la exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden	IPX-0	CAUTION: Protect from water / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Riskgefahr		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweis
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer	MEDIA	Media / Terreno di trasporto / Milieux / Medio / Medium
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	COL	Sample Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrenchen mit Schnitzverschluß	PRE REAG	Pretreatment Reagent / Reagente di Pretattamento / Réactif de prétraitement / Reactivo de prettamiento / Reagenz für die Vorbehandlung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampon de l'échantillon / Tampón de muestra / Probenpuffer	SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig		
TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß		

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.