alethia[®]

CMV DNA Amplification Assay

DNA Amplification Assay for the
Detection of Cytomegalovirus DNA

REF 481325

IVD

Rx Only

INTENDED USE

The Alethia CMV Assay Test System includes separately provided test kits for the Alethia CMV DNA Amplification Assay and the Alethia CMV External Control Reagents.

The Alethia CMV DNA Amplification Assay, performed on the Alethia instrument, is a qualitative, in vitro diagnostic test system for the direct detection of Cytomegalovirus (CMV) DNA in saliva samples from neonates younger than 21 days of age. The test is used as an aid in the diagnosis of congenital CMV infection. The results of this test should be used in conjunction with the results of other clinical findings.

Flocked swabs should be used to collect saliva from neonates. The swab can be collected dry, without viral transport media (VTM), or placed in no more than 1 mL VTM.

The Alethia CMV External Control Reagents are used as part of a routine quality control program to aid the user in detection of unexpected conditions that may lead to test errors. The external controls are intended for use with the Alethia CMV DNA Amplification Assay; the controls are not intended for use with other assays or systems.

SUMMARY AND EXPLANATION OF THE TEST

The Alethia CMV molecular assay is based on loop-mediated amplification (LAMP)^{1,2} technology. The assay targets a region of the Cytomegalovirus genome that is conserved across multiple CMV strains. The Alethia CMV target is a 194 base pair (bp) sequence of the Human herpesvirus 5 genome.

LAMP uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid solution. Reaction solution absorbance characteristics are monitored by the Alethia instrument. Changes in reaction solution turbidity created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia CMV kit includes Alethia CMV Test Devices and two separate CMV Buffers for sample processing and preparation. Saliva swab specimens are first treated with Buffer I (lysis buffer); viral capsids are disrupted and nucleic acids released. The lysate is then added to Buffer II (reaction buffer). The prepared sample contains nucleic acids for use with the Alethia CMV Test Device.

The Alethia CMV Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with Cytomegalovirus specific primers and a CONTROL chamber with human mitochondrial DNA-specific primers. Human mitochondrial DNA in saliva samples, and the human mitochondrial DNA-specific primers in the Test Device CONTROL chamber, function as the Internal Control for the assay. During specimen preparation, human mitochondrial DNA is liberated with the Cytomegalovirus DNA to allow for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors DNA amplification inhibition, assay reagent performance, and sample processing effectiveness. The Control target must be amplified and detected in the final reaction or the test is considered invalid and results are not reported.

The Alethia instrument monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{initial}, S_i) and at the assay Run End (Signal_{final}, S_f). The Alethia instrument calculates the change in light transmission between Run End and Run Start (S_f:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample assay results. TEST chamber S_f:S_i ratios less than 82% are reported as 'POSITIVE'. TEST chamber S_f:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. *Numerical values are not reported.*

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_f:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE or NEGATIVE). CONTROL chamber S_f:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. *Numerical values are not reported.*

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended, and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Cytomegalovirus (CMV) is a ubiquitous virus that infects people worldwide. Infections are generally acquired in youth with risk of infection increasing with age. Most infections are mild or unnoticed. CMV infection occurs through contact with body fluids from infected individuals. Infection in the USA adult population can range from approximately 42% to nearly 90%, depending on the demographic factor used in evaluation of antibody prevalence.³ In contrast to the mild infections in pediatric and adult patients, congenital CMV infections of newborns can be serious.⁴

CMV is the leading cause of congenital infections worldwide, and approximately 1-5% of births in developing countries are affected.⁴ Congenital CMV infection is a leading, nongenetic cause of sensorineural hearing loss.⁵ Identification of congenital CMV infection within the first three weeks of life may improve hearing and cognitive outcomes and is currently a focus of the American Academy of Pediatrics⁶ and the American Medical Association.⁷

Until recently, the gold standard method for diagnosing congenital CMV infection in many countries was CMV culture of urine or saliva collected within two weeks of birth.⁸ Studies indicate that salivary molecular methods, using dried or liquid saliva, offer excellent sensitivity and specificity for detecting congenital CMV.^{6,8-10} The Centers for Disease Control and Prevention lists saliva CMV tests as the standard method for diagnosing congenital CMV infection.¹¹

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Alethia **CMV Test Device (item #2835-01)**: Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) with Cytomegalovirus specific primers (TEST Chamber) and human mitochondrial DNA-specific primers (CONTROL Chamber).
2. Alethia **CMV Buffer I (item #2901-02)**: Lysis solution containing 0.2N sodium hydroxide and 1% Triton-X 100.
3. Alethia **CMV Buffer II (item #2779-02)**: Tris-buffer solution containing 0.09% azide as a preservative.

MATERIALS PROVIDED SEPARATELY

1. Alethia CMV External Control Kit; Meridian Bioscience, Inc. Catalog Number: 479880

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Specimen collection flocked swabs:
Dry Swab: Flocked swab in Dry Tube (e.g., Copan catalog 553C)

OR
Swab in Viral Transport Media: Flocked swab in Transport Media (Maximum Fill Volume: 1 mL; e.g., BD UVT catalog BD 220526 or Puritan UniTranz-RT UTM catalog UT116)

EQUIPMENT NOT PROVIDED

1. Interval Timer
2. Vortex mixer
3. Micropipette capable of dispensing 50 µL
4. Micropipette capable of dispensing 100 µL
5. Alethia Instrument; Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. CAUTION: The device contains material of human origin and should be handled as a potential carrier and transmitter of disease.
3. Follow Biosafety Level 2 and good laboratory practices during testing.¹² Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
5. Quality control programs for molecular testing laboratories, including workflow, and proper use and care of equipment, should be employed.¹³
6. Alethia CMV Test Devices contain lyophilized reagents. The protective pouches should not be opened until ready to perform the assay.
7. Alethia CMV Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used Alethia CMV Test Devices and tubes immediately after processing. Leave the Test Device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.
9. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343.3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
10. IMPORTANT: See SDS for additional safety and hazard information.

HAZARD AND PRECAUTIONARY STATEMENTS

 Buffer I	Signal Word Warning Hazard Statements H315 - Causes skin irritation H319 - Causes serious eye irritation Precautionary Statements - EU (§28, 1272/2008) P264 - Wash face, hands and any exposed skin thoroughly after handling. P280 - Wear eye protection/ face protection. P337 + P313 - If eye irritation persists: Get medical advice/attention. P321 – See SDS Section 4 or Section 11 for specific medical treatment information
---	--

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-30 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN COLLECTION AND PREPARATION

Sample type: Neonatal saliva swab. Saliva flocked swabs can be collected dry (without viral transport media) or placed in a viral transport media tube (no more than 1 mL).

Sample Collection: Collect saliva swab samples and place into a transport tube (with or without viral transport media) according to established laboratory methods. No special preparation of the neonate is required in order to collect the sample.

NOTE: Collect saliva swabs at least one hour after consumption of breastmilk.

Saliva swabs may be stored at room temperature (19-30 C), refrigerated (2-8 C), or frozen (≤ -20 C) after collection and during transportation to the laboratory. Samples should be tested as soon as possible but may be stored for up to 48 hours at 19-30 C, or 7 days refrigerated (2-8 C) prior to testing. Samples that will not be tested within 7 days should be frozen immediately at ≤ -20 C and may be stored up for to 14 days. Saliva samples on dry swabs or in VTM may be frozen and thawed up to 2 times after storage at ≤ -20 C prior to testing.

SPECIMEN PREPARATION

NOTE: Ensure the Alethia instrument is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia instrument Operator's Manual for further information regarding instrument set-up and operation.

NOTE: Ensure specimens are at room temperature (19-30 C) before specimen preparation.

NOTE: When vortexing solutions, care should be taken to avoid foaming to prevent splash, touch and aerosol contamination of target viral nucleic acids. Good laboratory practices should be adhered to when performing the Alethia CMV assay.

1. Swab preparation
 - **Dry Swab:** Remove cap from the Buffer I tube, and insert the swab head into the fluid. Rotate the swab head approximately 5 times to elute the saliva. Press the swab against the side of the tube to squeeze out fluid contents. Discard the swab and recap tube.
- OR**
 - **Swab in transport medium:** Remove the cap from the Buffer I tube. Using a micropipette, transfer 100 µL of transport medium to Buffer I tube and recap. Vortex tube for approximately 10 seconds to mix.
2. Incubate at 19-30 C for 2-3 minutes.
3. Remove the cap from the Buffer II tube. Using a micropipette, transfer 50 µL of Buffer I solution to the Buffer II tube. Recap the Buffer II tube.
4. Vortex the Buffer II tube for approximately 10 seconds to mix. Diluted samples may be held in Buffer II solution at 2-30 C for up to 3 hours. Samples must be warmed to room temperature (19-30 C) prior to preceding to TEST PROCEDURE.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single Alethia instrument run.

- For each sample, remove 1 Alethia CMV Test Device from its protective pouch. Carefully open the device by lifting the latch, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
- Using a micropipette, transfer 50 µL of the sample to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia CMV Test Device. Take care to not introduce air to the reaction mixtures. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
- Tap the device on the bench top or mix to remove air bubbles. Examine the Test Device for rehydration of the Control and Test Beads, and for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. No more than 10 minutes should elapse before initiating amplification and detection (i.e., amplification and detection should be initiated within 10 minutes of adding the sample to the Test Device).
- Insert the Alethia Test Device into the Alethia instrument and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Cytomegalovirus target DNA was detected
	NEGATIVE	No Cytomegalovirus DNA detected
	INVALID*	No reportable result. The test should be repeated: Samples can be retested from the Buffer II preparation within 3 hours of sample preparation.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use; Alethia instrument performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Run must be repeated. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia instrument performing correctly.
	INVALID	No reportable result. Run must be repeated. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia instrument well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat testing on the sample.

*Interpretation Notes

- Invalid results may occur as a result of inhibitory specimens, improper sample preparation, reagent failure, no human DNA in the specimen, instrument failure, or internal control failure.
- For VTM specimens, retesting may be performed with the original specimen if sufficient volume remains (see Specimen Collection and Preparation section for additional guidance).

QUALITY CONTROL

This test should be performed per applicable federal, state and local regulations or accrediting agencies.

- Each device contains an internal control that controls for amplification inhibition, assay reagents and sample processing effectiveness. Human mitochondrial DNA, which serves as the internal control DNA, is liberated from the saliva sample and processed through all steps of the procedure. Primers for amplification of human mitochondrial DNA are present in the Control Chamber of the Alethia Test Device.
- Good laboratory practice guidelines recommend the use of control materials.¹³ Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- Alethia CMV External Control Reagents are supplied separately (Catalog 479880). It is recommended that the reactivity of each new lot and each new shipment of Alethia CMV be verified on receipt or before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia CMV test kit should not be used in patient testing if the external controls do not produce the correct results.
- A separate Test Device must be used for each external control reagent.

EXPECTED VALUES

The observed prevalence of CMV from prospectively collected saliva specimens during this study was 0.3% (5/1,514). The results stratified by age and sex are provided below.

Sex/Age Group	Samples positive by CRM	Prevalence in Overall Prospectively Tested Population*
Female		
14-17 days	1	0.1%
Male		
0-3 days	4	0.3%
Total	5	0.3%

*Total number of prospectively tested specimens 1,514

LIMITATIONS OF THE PROCEDURE

- This product can only be used with the Alethia instrument.
- Samples containing mucin at concentrations >25 mg/mL may produce invalid results with the Alethia CMV assay.
- The results of the Alethia CMV assay are not intended to be used as the sole basis for diagnosis, treatment, or other patient management decisions.**
- The Alethia CMV assay is a qualitative assay and does not provide quantitative values or information about viral load.
- The detection of nucleic acids is dependent upon proper specimen collection, handling, transportation, storage, and preparation. Failure to observe proper procedure in any one of these steps may lead to incorrect results.
- Viral nucleic acid may persist *in vivo*, independent of viability. The Alethia CMV assay does not distinguish between viable and nonviable virus.
- As with all molecular based diagnostic tests: a) false negative results may occur from the presence of inhibitors, technical error, sample mix-up, or low numbers of organisms in the clinical specimen; b) false positive results may occur from the presence of cross-contamination by target virus nucleic acids or amplified product, and from non-specific signals.
- Saliva sample should be obtained at least one-hour after breast feeding to avoid contamination from cytomegalovirus which may result in a false positive result. Current treatment guidelines suggest waiting at least an hour after breastfeeding to obtain the sample.
- Viral transport media fill volumes greater than 1 mL are not acceptable for use with the Alethia CMV assay.
- Crossover contamination was observed for this device, it is recommended that good laboratory practices be adhered to when performing the Alethia CMV assay; gloves should be changed frequently to avoid carryover contamination.

SPECIFIC PERFORMANCE CHARACTERISTICS

CLINICAL PERFORMANCE

The Alethia CMV DNA Amplification Assay was evaluated from August 2017 to March 2018 at seven clinical study sites representing geographically distinct regions throughout the United States, Canada, Europe, and Australia. One-thousand five-hundred and fourteen (1,514) specimens were prospectively collected and tested with Alethia CMV assay (forty-seven of these specimens were frozen at ≤ -20 C after collection and tested later). The saliva swabs were collected at least one hour after breastfeeding.

For estimation of CRM positive percent agreement, thirty-four (34) archived specimens were also tested by Alethia CMV assay. The archived samples were de-identified samples previously evaluated from prospective clinical studies and found to have CMV infection. These samples were collected from infants less than 21 days of age and stored at -80 C after the completion of the initial testing. The status of breastfeeding time in relation to time of saliva samples collected was not available.

All samples were tested with the Alethia CMV assay at the study sites, then shipped to Meridian Bioscience, Inc. for composite reference method (CRM) testing. The CRM consisted of two manufacturer-developed and validated PCR assays. Samples positive by either PCR assay were further tested by bidirectional sequencing (BDS). Samples were considered positive when bi-directional sequencing results from either comparator PCR assay confirmed the presence of CMV DNA. Samples were considered negative when neither of the comparator PCR assays produced amplicon at the end of the 40-cycle amplification or BDS was negative.

Composite Reference Method

PCR #1	PCR #2	Bidirectional Sequencing	Composite Reference Method
+	+	+	+
+	+	-	-
+	-	+	+
+	-	-	-
-	+	+	+
-	+	-	-
-	-	n/a	-

Performance of the Alethia CMV assay was based on evaluation of 1,514 prospectively collected samples. From the prospectively collected samples; five samples were positive by CRM, 1,475 samples were negative by CRM (prevalence was 0.3% (5/1,514)), and 34 samples were invalid by the CRM and removed from the analysis (one of the samples was invalid by the Alethia CMV Assay). Among archived preselected positive samples, there were 34 samples positive by CRM.

Prospective Study

Alethia CMV	Composite Reference Method		Total
	Positive	Negative	
Positive	5	3	8
Negative	0	1472	1472
Total	5	1475	1480

	Estimate	95% CI
Positive Percent Agreement	100% (5/5)	56.7%; 100%
Negative Percent Agreement	99.8% (1,472/1,475)	99.4%; 99.9%

Preselected positive archived sample study:

34 archived preselected positive samples were evaluated in the study and confirmed positive by the CRM and the Alethia CMV Assay.

Alethia CMV	Composite Reference Method		Total
	Positive	Negative	
Positive	34	0	34
Negative	0	0	0
Total	34	0	34

	Estimate	95% CI
Positive Percent Agreement	100% (34/34)	89.9%; 100%
Negative Percent Agreement	N/A	N/A

Combined positive and negative percent agreement for the prospective study and the pre-selected positive sample study:

The positive samples from the prospective study (5) and archived pre-selected positive sample study (34) both demonstrated 100% agreement with the composite reference method and were thus presented combined in the table below.

	Estimate	95% CI
Positive Percent Agreement	100% (39/39)	91.0%; 100%
Negative Percent Agreement	99.8% (1,472/1,475)	99.4%; 99.9%

Invalid Results by Alethia CMV Assay

Twenty-seven (27) samples produced invalid results during initial testing with Alethia CMV assay. The rate of initial invalid results was 1.7% (27/1,548) and ranged from 0.7% to 4.1% at different clinical sites. After Alethia CMV assay re-testing, 26 samples had valid results. The rate of final invalid results was 0.06% (1/1,548) with 95% CI: 0.01%; 0.37%.

Additional Clinical Study Description:

Alethia CMV Assay Performance by Age (prospective and preselected positive sample types combined)

Age Group	Samples Positive by CRM			Samples Negative by CRM		
	Positive Percent Agreement	Alethia/ CRM	95% CI	Negative Percent Agreement	Alethia/CRM	95% CI
0-3 days	100%	15/15	79.6%; 100%	99.9%	1,351/1,353	99.5%; 100%
4-7 days	n/a	n/a	n/a	99.1%	108/109	95.0%; 99.8%
8-10 days	n/a	n/a	n/a	100%	4/4	51.0%; 100%
11-13 days	n/a	n/a	n/a	100%	2/2	34.2%; 100%
14-17 days	100%	1/1	20.7%; 100%	100%	6/6	61.0%; 100%
18-20 days	n/a	n/a	n/a	100%	1/1	20.7%; 100%
Known <21 days	100%	23/23	85.7%; 100%	n/a	n/a	n/a
Clinical Site Totals	100%	39/39	91.0%; 100%	99.8%	1,472/1,475	99.4%; 99.9%

Alethia CMV Assay Sex-Specific Performance (prospective and preselected positive sample types combined)

Sex	Samples Positive by CRM			Samples Negative by CRM		
	Positive Percent Agreement	Alethia/CRM	95% CI	Negative Percent Agreement	Alethia/CRM	95% CI
Female	100%	5/5	56.6%; 100%	99.7%	708/710	99.0%; 99.9%
Male	100%	11/11	74.1%; 100%	99.9%	763/764	99.3%; 100%
Unknown	100%	23/23	85.7%; 100%	100%	1/1	20.7%; 100%
Clinical Site Totals	100%	39/39	91.0%; 100%	99.8%	1,472/1,475	99.4%; 99.9%

ANALYTICAL SENSITIVITY

The analytical sensitivity, described as Limit of Detection (LoD), is the concentration at which the Alethia CMV assay has positive results at least 95% of the time. The LoD of the Alethia CMV assay was determined for Cytomegalovirus strain Merlin in a negative sample matrix. This negative matrix consisted of saliva samples collected from healthy CMV negative adults which were tested and confirmed negative by the Alethia CMV assay. Three kit lots of Alethia CMV assay reagents and eight Alethia instruments were used. A minimum of 6 dilutions with 20 replicates were tested for each lot. LoD was determined using Probit analysis. LoD was determined separately for the two sample types, dry swab and swab in VTM. LoD concentrations for each sample type are summarized below.

Sample Type	LoD (Copies/mL)
Dry Swab ¹	1,025
Swab in VTM ²	15,686

¹The CMV cp/mL concentration in saliva on the swab for dry swab samples is calculated by multiplying the CMV concentration in Buffer I by a factor of 4.75 (the 4.75-fold dilution occurs when 0.080 mL of saliva on the swab is added to 0.3 mL of Buffer I [380÷80 = 4.75]; 100% transfer is assumed).

²The CMV cp/mL concentration in saliva on the swab for VTM samples is calculated by multiplying the CMV concentration in VTM by a factor of 13.5 (the 13.5-fold dilution occurs when 0.080 mL of saliva on the swab is placed in 1 mL of VTM [1080÷80 = 13.5]; 100% transfer is assumed).

ASSAY INCLUSIVITY

Samples containing 3 additional CMV strains other than CMV strain Merlin were evaluated with the Alethia CMV assay. These strains include Toledo, Towne, and AD-169. Quantified strains were diluted in simulated negative clinical matrix to approximately 2-3X LoD for both VTM and dry swab and tested in triplicate. All strains tested for both sample types produced positive results with Alethia CMV.

CMV Strain Tested	Dry Swab Samples		VTM Swab Samples	
	Test concentration (copies/mL)	Test Results (n _{pos} /n _{total})	Test concentration (copies/mL)	Test Results (n _{pos} /n _{total})
AD-169	2,493	3/3	45,077	3/3
Toledo	2,493	3/3	45,077	3/3
Towne	2,493	3/3	45,077	3/3

REPRODUCIBILITY

Panels of two sample types, dry swab and VTM, were supplied to three laboratories for this reproducibility study. The panels included contrived CMV samples manufactured as moderate-positive, low-positive samples, and high-negative samples (30 replicates per site for each sample). The panel also included one true negative sample (10 replicates per site). Positive and Negative Controls were tested with each panel also (10 replicates per site). Testing was performed by different operators at each site on the same day for five days. Three lots of Alethia CMV kits and six Alethia instruments were used in this study. Mean values, repeatability, between-operator, between-day, and between-site components of variance, percent of positive and negative results are provided in the tables below for each sample type tested in the reproducibility study.

Reproducibility for Dry Swab Samples

Saliva Samples on Dry Swabs														
	N	Mean	Repeatability		Between-operators ¹		Between-day		Between-site ²		Reproducibility		%Pos results	%Neg results
			SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV		
Moderate Positive	90	62.26	2.93	4.7%	1.69	2.7%	0.00	0.0%	0.91	1.5%	3.50	5.6%	100%	0%
Low Positive	90	62.38	4.44	7.1%	1.96	3.1%	1.11	1.8%	0.00	0.0%	4.98	8.0%	98.9%	1.1%
High Negative	90	99.05	6.68	6.7%	0.97	1.0%	0.00	0.0%	0.23	0.2%	6.76	6.8%	3.3%	96.7%
True Negative ³	30	100.13	0.69	0.7%	n/a	n/a	0.31	0.3%	0.40	0.4%	0.86	0.9%	0%	100%
Positive Control ³	30	61.84	2.72	4.4%	n/a	n/a	0.00	0.0%	1.11	1.8%	2.94	4.8%	100%	0%
Negative Control ³	30	100.23	2.06	2.1%	n/a	n/a	0.17	0.2%	0.22	0.2%	2.08	2.1%	0%	100%

¹ Includes between-operator and between-instrument components

² Includes between-site and between-kit lot components

³ Samples were run at each site for 5 days with 2 runs per day and 1 replicate per run.

Reproducibility for VTM samples

Saliva Samples in VTM														
			Repeatability		Between-operators ¹		Between-day		Between-site ²		Reproducibility		%Pos results	%Neg results
	N	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV		
Moderate Positive	90	61.99	3.18	5.1%	1.02	1.6%	0.00	0.0%	1.20	1.9%	3.55	5.7%	100%	0%
Low Positive	90	61.71	2.37	3.8%	1.43	2.3%	0.00	0.0%	0.67	1.1%	2.85	4.6%	100%	0%
High Negative	90	100.11	4.82	4.8%	0.00	0.0%	0.74	0.7%	0.00	0.0%	4.88	4.9%	1.1%	98.9%
True Negative ³	30	99.20	7.71	7.8%	N/A	N/A	0.00	0.0%	0.00	0.0%	7.71	7.8%	3.3%	96.7%
Positive Control ³	30	62.22	2.39	3.8%	N/A	N/A	0.00	0.0%	0.49	0.8%	2.44	3.9%	100%	0%
Negative Control ³	30	100.30	2.02	2.0%	N/A	N/A	0.00	0.0%	0.12	0.1%	2.02	2.0%	0%	100%

¹ Includes between-operator and between-instrument components

² Includes between-site and between-kit lot components

³ Samples were run at each site for 5 days with 2 runs per day and 1 replicate per run.

WITHIN-LABORATORY PRECISION STUDY

Panels of two sample types, dry swab and VTM, were tested at one site (internal) over 6 days. The panels included contrived CMV samples manufactured as moderate-positive, low-positive (around 2X LoD), and high-negative samples. The panel also included one true negative sample, Positive and Negative Controls. Three kit lots were used during the study, one lot per day of testing. Each kit lot was tested twice over the 6-day testing period. Mean values, repeatability, between-operator, between-day, and between-kit lot components of variance, percent of positive and negative results are provided in the tables below for each sample type tested in the precision study.

Within-Laboratory Precision for Dry Swab Samples

Saliva Samples on Dry Swabs														
			Repeatability ¹		Between-operators ¹		Between-day		Between-kit lot		Reproducibility		%Pos results	%Neg results
	N	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV		
Moderate Positive	48	60.29	1.99	3.3%	0.00	0.0%	0.00	0.0%	0.74	1.2%	2.12	3.5%	100%	0%
Low Positive	48	60.23	1.70	2.8%	0.62	1.0%	0.00	0.0%	0.74	1.2%	1.95	3.2%	100%	0%
High Negative	48	98.26	6.46	6.6%	4.36	4.4%	1.42	1.4%	0.80	0.8%	7.96	8.1%	4.2%	95.8%
True Negative	48	99.94	0.61	0.6%	0.00	0.0%	0.12	0.1%	0.14	0.1%	0.64	0.6%	0%	100%
Positive Control	24	60.57	1.96	3.2%	0.00	0.0%	1.47	2.4%	0.00	0.0%	2.45	4.0%	100%	0%
Negative Control	24	99.73	0.29	0.3%	0.00	0.0%	0.18	0.2%	0.16	0.2%	0.38	0.4%	0%	100%

¹ Includes between-operator and between-instrument components

Within-Laboratory Precision for VTM Samples

Saliva Samples in VTM														
			Repeatability		Between-operators ¹		Between-day		Between-kit lot		Reproducibility		%Pos results	%Neg results
	N	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV		
Moderate Positive	48	60.62	1.36	2.2%	0.91	1.5%	0.00	0.0%	0.77	1.3%	1.81	3.0%	100%	0%
Low Positive	48	60.55	1.67	2.8%	0.89	1.5%	0.72	1.2%	0.00	0.0%	2.02	3.3%	100%	0%
High Negative	48	99.91	0.64	0.6%	0.00	0.0%	0.22	0.2%	0.00	0.0%	0.68	0.7%	0%	100%
True Negative	48	100.00	0.84	0.8%	0.28	0.3%	0.00	0.0%	0.00	0.0%	0.89	0.9%	0%	100%
Positive Control	24	60.28	2.12	3.5%	0.00	0.0%	1.09	1.8%	0.00	0.0%	2.38	4.0%	100%	0%
Negative Control	24	100.29	1.53	1.5%	0.28	0.3%	0.00	0.0%	0.00	0.0%	1.56	1.6%	0%	100%

¹ Includes between-operator and between-instrument components

CROSS REACTIVITY STUDIES

Cross-reactivity studies were performed with a panel of 40 microorganisms and human genomic DNA, each diluted in dry swab simulated negative clinical matrix. Microorganisms selected include those with genetic similarity to CMV and those likely to be present in the oral cavity of neonates. Microorganisms were diluted in dry swab simulated negative clinical matrix to the indicated concentration (see table below) and tested in triplicate. No cross-reactivity with the Alethia CMV assay was observed.

Microorganism	Test concentration	Test Results (n _{pos} /n _{total})	Microorganism	Test concentration	Test Results (n _{pos} /n _{total})
<i>Acinetobacter baumannii</i>	1.2x10 ⁷ CFU/mL	0/3	Adenovirus	3.80x10 ⁵ TCID ₅₀ /mL	0/3
<i>Actinomyces odontolyticus</i>	1.2x10 ⁷ CFU/mL	0/3	Coronavirus	2.19x10 ⁵ TCID ₅₀ /mL	0/3
<i>Bordetella pertussis</i>	1.2x10 ⁷ CFU/mL	0/3	Coxsackievirus	4.07x10 ⁶ TCID ₅₀ /mL	0/3
<i>Candida albicans</i>	1.2x10 ⁷ CFU/mL	0/3	Enterovirus 71	1.26x10 ⁵ TCID ₅₀ /mL	0/3
<i>Escherichia coli</i> ATCC 35218	1.2x10 ⁷ CFU/mL	0/3	Epstein Barr Virus	3.39x10 ⁸ cp/mL	0/3
<i>Fusobacterium nucleatum</i>	1.2x10 ⁷ CFU/mL	0/3	Herpes Simplex Virus 1	9.5x10 ⁵ TCID ₅₀ /mL	0/3
<i>Haemophilus influenzae</i>	1.2x10 ⁷ CFU/mL	0/3	Herpes Simplex Virus 2	1.3x10 ⁵ TCID ₅₀ /mL	0/3
<i>Haemophilus parainfluenzae</i>	1.2x10 ⁷ CFU/mL	0/3	Human herpesvirus 6B	6.16x10 ⁷ cp/mL	0/3
<i>Moraxella catarrhalis</i>	1.2x10 ⁷ CFU/mL	0/3	Human herpesvirus 7	3.80x10 ⁵ TCID ₅₀ /mL	0/3
<i>Mycoplasma pneumoniae</i>	3.70x10 ⁷ CCU/mL	0/3	Human herpesvirus 8	2.13x10 ⁸ cp/mL	0/3
<i>Porphyromonas gingivalis</i>	1.2x10 ⁷ CFU/mL	0/3	Human metapneumovirus	6.61x10 ⁵ TCID ₅₀ /mL	0/3
<i>Pseudomonas aeruginosa</i>	1.2x10 ⁷ CFU/mL	0/3	Influenza A	3.80x10 ⁵ TCID ₅₀ /mL	0/3
<i>Staphylococcus aureus</i>	1.2x10 ⁷ CFU/mL	0/3	Influenza B	4.57x10 ⁵ TCID ₅₀ /mL	0/3
<i>Staphylococcus epidermidis</i>	1.2x10 ⁷ CFU/mL	0/3	Parainfluenza virus 1	1.95x10 ⁶ TCID ₅₀ /mL	0/3
<i>Streptococcus agalactiae</i> (GBS)	1.2x10 ⁷ CFU/mL	0/3	Parainfluenza virus 2	5.89x10 ⁶ TCID ₅₀ /mL	0/3
<i>Streptococcus anginosus</i> (Group F)	1.2x10 ⁷ CFU/mL	0/3	Parainfluenza virus 3	2.19x10 ⁵ TCID ₅₀ /mL	0/3
<i>Streptococcus mitis</i>	1.2x10 ⁷ CFU/mL	0/3	Respiratory syncytial virus A	3.2x10 ⁵ TCID ₅₀ /mL	0/3
<i>Streptococcus oralis</i>	1.2x10 ⁷ CFU/mL	0/3	Respiratory syncytial virus B	4.6x10 ⁵ TCID ₅₀ /mL	0/3
<i>Streptococcus salivarius</i>	1.2x10 ⁷ CFU/mL	0/3	Rhinovirus	1.51x10 ⁵ TCID ₅₀ /mL	0/3
<i>Streptococcus sanguinis</i>	1.2x10 ⁷ CFU/mL	0/3	Varicella Zoster Virus	3.36x10 ⁹ cp/mL	0/3
Human Genomic DNA	6.18x10 ⁶ cp/mL	0/3	n/a	n/a	n/a

TESTS FOR MICROBIAL INTERFERENCE

Microbial interference studies were performed with a panel of 40 microorganisms and human genomic DNA, each diluted in dry swab simulated positive clinical matrix. Each sample tested contained CMV at a concentration of 3X LoD and microorganism or human genomic DNA at the test concentration indicated above (see table in Cross-reactivity section for microorganisms tested and final test concentrations). Once prepared, each sample was tested in triplicate. No microbial interference with the Alethia CMV assay was observed (i.e., all microorganisms diluted in dry swab simulated positive clinical matrix produced 3/3 positive test results with the Alethia CMV assay).

TESTS FOR INTERFERING SUBSTANCES

Interference testing was performed in the presence of chemical and biological substances introduced directly into contrived CMV low positive (3X LoD) and negative samples. Two unique positive samples and one negative sample were tested in triplicate. The substances tested, the concentrations evaluated, and test results are shown in the following table. No interference was observed with the following substances with the Alethia CMV assay at the specified test concentrations (i.e., all positive replicates tested produced positive Alethia CMV test results and all negative replicates tested produced negative Alethia CMV test results).

Substance Tested	Test concentration	Test Results (n _{detected} /n _{total})		
		CMV Negative Sample	CMV Positive Sample #1	CMV Positive Sample #2
Infants' Pain & Fever (Acetaminophen)	0.2 mg/mL	0/3	3/3	3/3
Acetylsalicylic acid	0.65 mg/mL	0/3	3/3	3/3
Caffeine	0.06 mg/mL	0/3	3/3	3/3
Enfamil™ Fer-In-Sol® (Ferrous Sulfate)	1.5 mg/mL	0/3	3/3	3/3
Enfamil Premium® Infant Formula Newborn	10% v/v	0/3	3/3	3/3
Infants' Mylicon® Gas Relief (Simethicone)	2 mg/0.3 mL	0/3	3/3	3/3
Gaviscon® infant (Sodium alginate)	1.2 mg/mL	0/3	3/3	3/3
Magnesium alginate	0.467 mg/mL	0/3	3/3	3/3
Infants' Ibuprofen	0.5 mg/mL	0/3	3/3	3/3
Enfamil™ Poly-Vi-Sol®	8% v/v	0/3	3/3	3/3
Little Remedies® Saline spray/drops	10% v/v	0/3	3/3	3/3
Methadone	0.002 mg/mL	0/3	3/3	3/3
Morphine sulphate	0.0005 mg/mL	0/3	3/3	3/3
Nystatin	1000 U/mL	0/3	3/3	3/3
Prednisone	0.0003 mg/mL	0/3	3/3	3/3
Casein	10 mg/mL	0/3	3/3	3/3
Mucin (from porcine stomach – Type II)*	25 mg/mL	0/3	3/3	3/3
White blood cells	10% v/v	0/3	3/3	3/3
Whole blood	10% v/v	0/3	3/3	3/3

*When mucin was tested at a concentration of 50 mg/mL, the negative sample tested negative in 3/3 replicates and CMV positive sample #1 tested positive in 3/3 replicates. CMV positive sample #2 produced two invalid test results, and one positive test result for the three replicates tested. The mucin concentration in the samples was reduced to 25 mg/mL and testing was repeated, producing the results shown in the table. No interference was observed with a mucin concentration of 25 mg/mL.

CARRYOVER CONTAMINATION STUDIES

A series of alternating true negative and high viral load samples (30X LoD) were processed using the Alethia CMV assay. A total of 30 high positive and 30 true negative samples were processed. Six runs were conducted with 10 samples each (5 high positive and 5 true negative samples). All samples tested produced expected results with all positive samples producing positive Alethia CMV test results, and all negative samples producing negative Alethia CMV test results. No carryover was observed. Workflow studies revealed that foam present after vortexing samples presents an aerosol contamination risk (i.e., when bubbles break), and/or a transfer contamination risk if gloves come in contact with liquids. It is recommended that good laboratory practices be adhered to when performing the Alethia CMV assay; gloves should be changed frequently.

SAMPLE STABILITY STUDY SUMMARY

A panel of contrived CMV low positive (3X LoD) and negative samples of both types (dry swab and VTM) were prepared and stored in aliquots as indicated below. Aliquots of each sample were removed from storage at specified intervals and allowed to equilibrate to room temperature prior to testing in triplicate. For the storage parameters indicated, all positive samples produced positive Alethia CMV results, and all negative samples produced negative Alethia CMV results.

Storage Condition	Acceptable Storage Parameters
Room Temperature (19-30 C)	48 hours
Refrigerated (2-8 C)	7 days
Frozen (≤ -20 C)	14 days, may be frozen and thawed up to 2 times

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REFERENCES

1. Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002;16:223-29.
2. Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. *J Biochem Biophys* 2004;59:145-47.
3. Bate SL et al. Cytomegalovirus seroprevalence in the United States: The National Health and Nutrition Examination Surveys, 1988-2004. *CID* 2010;50:1439-1447.
4. Manicklal S et al. The silent global burden of congenital cytomegalovirus. *Clin Microbiol Rev* 2013;26:86-102.
5. Boppana S et al. Saliva polymerase-chain-reaction assay for Cytomegalovirus screening in newborns. *NEJM* 2011;364:2111-2118.
6. American Academy of Pediatricians, Congenital Cytomegalovirus (CMV) infection and hearing loss. www.aap.org/aap-health-initiatives/PEHDIC/Documents.pdf.
7. Demmler-Harrison GJ. Editorial: Congenital cytomegalovirus infection. The elephant in our living room. jamanetwork.com/journals/jamapediatrics/fullarticle/2557386.
8. Ross SA et al. Detection of congenital cytomegalovirus infection by real-time polymerase chain reaction analysis of saliva or urine specimens. *JID* 2014;210:1415-1418.
9. Yamamoto Y et al. Is saliva as reliable as urine for detection of cytomegalovirus DNA for neonatal screening of congenital CMV infection? *J Clin Virol* 2006;36:228-230.
10. Bélec L, Brogan TV. Real-time PCR-based testing of saliva for cytomegalovirus at birth. *Expert Rev anti-Infect Ther* 2011;9:1119-1124.
11. Centers for Disease Control: Laboratory diagnosis of CMV infection for persons \geq 12 months of age. <http://cdc.gov/cmvc/clinical/lab-tests.html>.
12. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
13. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

SN11118


REV. 03/26

 <p>Manufactured By</p>	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
<p>EC REP</p>	<p>Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders:</p> <ul style="list-style-type: none"> • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu
<p>UK Authorised Representative</p>	<p>Launch Diagnostics Limited Lakeview West Crossways Business Park Galleon Boulevard Dartford Kent DA2 6QE United Kingdom</p>
<p>CH REP</p>	<p>MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland</p>

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	Meridian products carrying the European Conformity (CE) mark fulfil the requirements of Directive 98/79/EC or the Regulation 2017/746 on in-vitro diagnostic medical devices / I prodotti Meridian recanti il marchio di Conformità Europea (CE) soddisfano i requisiti della Direttiva 98/79/CE o del Regolamento 2017/746 sui dispositivi medico-diagnostici in vitro / Les produits Meridian portant la marque de Conformité européenne (CE) sont conformes aux exigences de la Directive 98/79/CE ou du Règlement 2017/746 relatifs aux dispositifs médicaux de diagnostic in vitro / Los productos de Meridian que llevan la marca de conformidad europea (CE) cumplen los requisitos de la Directiva 98/79/CE o del Reglamento 2017/746 sobre los productos sanitarios para diagnóstico in vitro / Produkte von Meridian mit der CE-Kennzeichnung erfüllen die Anforderungen der EU-Richtlinie 98/79/EG bzw. der Verordnung (EU) 2017/746 über In-vitro-Diagnostika	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
		STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Stérilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Caution, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Stérilisation par oxyde d'éthylène / Esterilizado por oxidode etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femminile / De sexe féminin / Hembra / Frau		Recycle – do not dispose of as general waste / Riciclare – non eliminare come rifiuto generico / Recycler – ne pas jeter dans une poubelle / Recycle – no deseché como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung	HT TUBE	Heat Treatment Tube / Provetta per il Trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrrchen zur Hitzebearbeitung
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Eviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden		HOT SURFACE: Keep hands Away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	IPX-0	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	MEDIA	Media / Terreno di trasporto / Milieux / Medio / Medium
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrrchen mit Schnappverschluss	COL	Sample Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampón de l'échantillon / Tampón de muestra / Probenpuffer	PRE REAG	Pretreatment Reagent / Reagente di Pretrattamento / Réactif de prétraitement / Reactivo de pretratamiento / Reagenz für die Vorbehandlung
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung
TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß	CH REP	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.