
ImmunoCard® STAT₁ Strep A Test

ImmunoCard[®] STAT! Strep A Test

CLIA Complexity:
Waived

REF 755250

IVD In Vitro Diagnostic Medical Device

Rx Only

FOR LABORATORY AND PROFESSIONAL IN VITRO DIAGNOSTIC USE ONLY.

INTENDED USE

The ImmunoCard STAT! Strep A Test is intended for the qualitative detection of Group A Streptococcal antigen from throat swabs or confirmation of presumptive Group A Streptococcal colonies recovered from culture.

SUMMARY AND EXPLANATION OF TEST

Group A Streptococcus is one of the most important causes of acute upper respiratory tract infection. Early diagnosis and treatment of Group A Streptococcal pharyngitis has been shown to reduce the severity of symptoms and further complications such as rheumatic fever and glomerulonephritis.¹ Conventional identification procedures for Group A Streptococcus from throat swabs involve the isolation and subsequent identification of viable pathogens by techniques that require 24 to 48 hours or longer.² The ImmunoCard STAT! Strep A Test detects either viable or nonviable organisms directly from a throat swab, providing results within five minutes.

BIOLOGICAL PRINCIPLES

The ImmunoCard STAT! Strep A Test uses color immunochromatographic dipstick technology with rabbit antibodies coated on the nitrocellulose membrane. In the test procedure, a throat swab is subjected to a chemical extraction of a carbohydrate antigen unique to Group A Streptococcus. The Test Stick is then placed in the extraction mixture and the mixture migrates along the membrane. If Group A Streptococcus is present in the sample, it will form a complex with the anti-Group A Streptococcus antibody conjugated color particles. The complex will then be bound by the anti-Group A Streptococcus capture antibody and a visible blue Test Line will appear to indicate a positive result.

REAGENTS/MATERIALS PROVIDED AND STORAGE

- 50 Test Sticks
- 50 Test Tubes
- 50 Sterile Swabs
- 1 Reagent 1 (2 M Sodium Nitrite)
- 1 Reagent 2 (0.3 M Acetic Acid)
- 1 Positive Control (Nonviable Group A Streptococci, 0.1% Sodium Azide)
- 1 Negative Control (Nonviable Group C Streptococci, 0.1% Sodium Azide)
- 1 Directional Insert

Note: Extra components (swabs, tubes) have been provided for your convenience

Store Test Sticks and reagents tightly capped at 15-30 C.
Do not use Test Sticks or reagents after expiration date.

MATERIALS REQUIRED BUT NOT PROVIDED

A timer or watch.

PRECAUTIONS

For in vitro diagnostic use.

Follow your laboratory safety guidelines in the collection, handling, storage and disposal of controls, patient specimens and all items exposed to patient specimens.³

Reagent 2 contains an acid. If the solution comes in contact with the skin or eyes, flush with large volumes of water.

The Positive and Negative Controls contain 0.1% sodium azide which may react with lead or copper plumbing to form potentially explosive metal azide. For sites permitted to dispose of material down a sink: large quantities of water must be used to flush discarded control material down a sink. Do not interchange or mix components from different kit lots.

HAZARDS AND PRECAUTIONARY STATEMENTS

Component (s)	Pictogram	Signal Word	Hazardous Ingredients
ImmunoCard STAT! [®] Strep A REAG 1		Warning	sodium nitrite (CAS No) 7632-00-0
Hazard statements	H302 - Harmful if swallowed		
Precautionary statements	P264 - Wash hands, forearms and face thoroughly after handling. P270 - Do not eat, drink or smoke when using this product. P301+P312 - If swallowed: Call a poison center or doctor if you feel unwell. P330 - Rinse mouth. P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.		
ImmunoCard STAT! [®] Strep A REAG 2		Warning	Acetic acid (CAS No) 64-19-7
Hazard statements	H315 - Causes skin irritation H319 - Causes serious eye irritation		
Precautionary statements	P264 - Wash hands, forearms and face thoroughly after handling. P280 - Wear protective gloves/protective clothing/eye protection/face protection. P302+P352 - If on skin: Wash with plenty of water. P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P321 - Specific treatment (see supplemental first aid instruction on this label). P332+P313 - If skin irritation occurs: Get medical advice/attention. P337+P313 - If eye irritation persists: Get medical advice/attention. P362+P364 - Take off contaminated clothing and wash it before reuse.		

SPECIMEN COLLECTION AND PREPARATION

Collect specimens with a sterile swab from the tonsils and/or back of the throat² taking care to avoid the teeth, gums, tongue or cheek surfaces.

- Do not use swabs with cotton tips, wooden shafts or calcium alginate swabs.
- Do not use a collection system that contains charcoal or semisolid transport media.
- If your lab requires a culture result as well as the ImmunoCard STAT! Strep A Test result, streak the culture plate with the swab before starting the ImmunoCard STAT! Strep A Test procedure as the extraction reagents will cause the specimen to become nonviable.
- Process the swab as soon as possible after collecting the specimen. If you do not perform the ImmunoCard STAT! Strep A test immediately, store the swabs either at room temperature or refrigerated for up to 48 hours. The swabs and the test kit must be at room temperature prior to running the test.

Sample Transport:

- Because the performance characteristics of this product were established with the sterile rayon swabs supplied with the kit, we recommend using these swabs to assure optimal performance. In the U.S., you may purchase the kit swabs in a double swab/dry tube format as an accessory (Contact Meridian Bioscience Technical Services for more information).
- Because the test does not require live organisms for processing, a rayon transport swab containing Stuart's or Amies media may also be used; however, swabs from other suppliers have not been validated.


CULTURE CONFIRMATION

The ImmunoCard STAT! Strep A Test can also be used to confirm the identification of Group A Streptococcus on blood agar plates. The plates must be less than 72 hours old. Lightly touch 1-3 suspect colonies (showing characteristic beta hemolysis) using a sterile swab. Do not sweep the plate. Follow the instructions in the TEST PROCEDURE section to test the swab.


TEST PROCEDURE


- Just before testing, add 3 drops Reagent 1 (pink) and 3 drops Reagent 2 to the Test Tube (the solution should turn light yellow).


- Immediately put the swab into the Tube.
- Vigorously mix the solution by rotating the swab forcefully against the side of the Tube at least ten (10) times. Best results are obtained when the specimen is vigorously extracted in the solution.
- Let stand for 1 minute.

- Express as much liquid as possible from the swab **by squeezing the sides of the tube as the swab is withdrawn.**
- Discard the swab.

- Remove Test Stick(s) from the container; re-cap container immediately.
- Place the Absorbent End of the Test Stick into the extracted sample.
- Read results at 5 minutes. Positive results may be read as soon as the red Control Line appears.
- Results are invalid after the stated read time. The use of a timer is recommended.**

INTERPRETATION OF TEST RESULTS

Notes

A blue or red line which appears uneven in color density is considered a valid result. In cases of moderate or high positive specimens, some blue color behind the Test Line may be seen; as long as the Test Line and Control Line are visible, the results are valid.

Positive

A blue Test Line and a red Control Line is a positive result for the detection of Group A Streptococcus antigen. Note that the blue line can be any shade of blue and can be lighter or darker than the line in the picture.

Negative

A red Control Line but no blue Test Line is a presumptive negative result.

Invalid

If no red Control Line appears or background color makes reading the red Control Line impossible, the result is invalid. If this occurs, repeat the test on a new Test Stick or contact Meridian Bioscience Technical Services.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

Internal Procedural Controls

The ImmunoCard STAT! Strep A Test provides three levels of procedural controls with each test run:

- The color of the liquid changes from pink to light yellow as you add Extraction Reagent 2 to Extraction Reagent 1. This is an internal extraction reagent control. The color change means that you mixed the extraction reagents properly. The color change also means that the reagents are functioning properly.
- The red Control Line is an internal positive procedural control. The Test Stick must absorb the proper amount of sample and the Test Stick must be working properly for the red Control Line to appear. For the Test Stick to be working properly, the capillary flow must occur.
- A clear background is an internal background negative procedural control. If no interfering substances are in the specimen and the Test Stick is working properly, the background in the Control Line area will clear. A discernible result will be seen.

If the red Control Line does not appear, the test may be invalid. If the background does not clear and interferes with the test result, the test may be invalid. If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

External Quality Control Testing

Each kit contains Positive and Negative Control material. The Controls are for external quality control testing. Use the Controls to test that the extraction reagents and the Test Sticks are working. Also use the Controls to test that you are able to correctly perform the test procedure. If you choose, you may use Group A and non Group A Streptococcus ATCC reference strains as controls. Some commercial

controls may contain interfering additives. Therefore Meridian Bioscience recommends that you do not use other commercial controls with the ImmunoCard STAT! Strep A Test.

Quality Control requirements should be established in accordance with local, state and federal regulations or accreditation requirements. Minimally, Meridian Bioscience recommends that positive and negative external controls be run with each new lot and with each new untrained operator.

QC Testing Procedure:

- Dispense 3 drops Reagent 1 and 3 drops Reagent 2 into Test Tube.
- Vigorously mix the control contents. Add 1 free falling drop of Control from dropper bottle.
- Place a clean swab into the Tube.
- Continue as you would for a patient sample, as instructed in the PROCEDURE section.

EXPECTED VALUES

Approximately 19% of all upper respiratory tract infections are caused by Group A Streptococci.⁵ Streptococcal pharyngitis displays a seasonal variation and is most prevalent during winter and early spring. The highest incidence of this disease is found in crowded populations such as military bases and in school-age children.⁶

LIMITATIONS OF THE PROCEDURE

- The ImmunoCard STAT! Strep A Test has been categorized as CLIA waived only for the application of qualitative detection of Group A Streptococcal Antigen from throat swabs. The application for the confirmation of presumptive Group A Streptococcal colonies recovered from culture is not waived.
- The results obtained with this kit yield data that must be used only as an adjunct to other information available to the physician. The ImmunoCard STAT! Strep A Test is a qualitative test for the detection of Group A Streptococcal antigen. This test does not differentiate between viable and nonviable Group A Streptococci.
- The ImmunoCard STAT! Strep A Test should be used only with throat swabs or colonies taken directly from a plate. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper throat swab specimens must be obtained.
- This test does not differentiate between carriers and acute infection. Pharyngitis may be caused by organisms other than Group A Streptococcus.^{1,2}
- A negative result may be obtained if the specimen is inadequate or antigen concentration is below the sensitivity of the test.
- The American Academy of Pediatrics states⁴: "Several rapid diagnostic tests for GAS pharyngitis are available. The specificities of these tests generally are very high, but the reported sensitivities vary considerably. As with throat cultures, the accuracy of these tests is most dependent on the quality of the throat swab specimen, which must contain pharyngeal and tonsillar secretions, and on the experience of the person who is performing the test. Therefore, when a patient suspected of having GAS pharyngitis has a negative rapid streptococcal test, a throat culture should be obtained to ensure that the patient does not have GAS infection." It also states: "Cultures that are negative for GAS infection after 24 hours should be incubated for a second day to optimize isolation of GAS."

SPECIFIC PERFORMANCE CHARACTERISTICS

In a multi-center evaluation, a total of 639 throat swabs were collected from patients presenting with pharyngitis. Each swab was inoculated to a sheep blood agar plate, then tested by the ImmunoCard STAT! Strep A Test. Plates were incubated for 18-24 hours at 35-37 °C at 5-10% CO₂ with a Bacitracin disk. Presumptive GAS colonies were confirmed with commercially available Strep A testing kits.

Of the 639 total specimens, 464 were found to be negative by culture and 454 were also negative by the ImmunoCard STAT! Strep A Test, for a specificity of 97.8%. Of the 175 specimens found to be positive by culture, 168 were also positive by the ImmunoCard STAT! Strep A Test, for a sensitivity of 96.0%. The 95% confidence intervals were calculated to be 96.6-99.0% for specificity and 94.4-97.6% for sensitivity. Overall agreement between culture and the ImmunoCard STAT! Strep A Test was 97.3% (622/639). The results are summarized below:

Culture Classification	ImmunoCard STAT! / Culture	% Correct
Negative (Specificity)	454/464	97.8%
1+ (≤ 10 colonies)	3/6	50.0%
2+ (11-50 colonies)	9/13	69.2%
3+ (> 50 colonies)	44/44	100%
4+ (predominant growth)	112/112	100%
Total Positive (Sensitivity)	168/175	96.0%
Total (Overall Agreement)	622/639	97.3%

In addition, the ImmunoCard STAT! Strep A Test was used to confirm the identification of Group A Streptococcus on blood agar plates. As a culture confirmation test, the ImmunoCard STAT! Strep A Test was 100% sensitive (62/62) and 100% specific (39/39).

Crossreactivity

The following organisms tested at levels of approximately 1×10^8 organisms/test were all found to be negative when tested with the ImmunoCard STAT! Strep A Test.

<i>Streptococcus</i> Group B	<i>Staphylococcus aureus</i>	<i>Bordetella pertussis</i>
<i>Streptococcus</i> Group C	<i>Staphylococcus epidermidis</i>	<i>Neisseria meningitidis</i>
<i>Streptococcus</i> Group F	<i>Corynebacterium diphtheriae</i>	<i>Neisseria gonorrhoeae</i>
<i>Streptococcus</i> Group G	<i>Serratia marcescens</i>	<i>Neisseria sicca</i>
<i>Streptococcus pneumoniae</i>	<i>Candida albicans</i>	<i>Neisseria subflava</i>
<i>Streptococcus sanguis</i>	<i>Klebsiella pneumoniae</i>	<i>Branhamella catarrhalis</i>
<i>Streptococcus mutans</i>	<i>Pseudomonas aeruginosa</i>	<i>Hemophilus influenzae</i>
<i>Enterococcus faecalis</i>		

POL Studies

An evaluation of the ImmunoCard STAT! Strep A Test was conducted at three physicians offices where testing was performed by personnel with diverse educational backgrounds. Each site tested the randomly coded panel consisting of negative (6), low positive (3) and moderate positive (3) specimens for three days. The results obtained had >99% agreement (107/108) with the expected results.

RE-ORDER

No.755250 ImmunoCard STAT! Strep A Test (50 Tests)

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Nr. 755250 ImmunoCard STAT! Strep A Test (50 Tests)

REFERENCES

1. Youmans GP, Paterson PY, and Sommers HM. Upper respiratory tract infection: general considerations, in the biologic and clinical basis of infectious diseases, W. B. Saunders Co., Philadelphia, 1980; 177-183.
2. Balows A, Hausler WJ, Herrmann KL, Isenberg HD, Shadomy HJ, Eds., Faklam, RR and Washington JA, Streptococcus and related catalase-negative gram-positive cocci, in Manual of Clinical Microbiology, 5th Edition, Am. Society of Microbiology, Washington, D.C., 1991; 238-257.
3. CDC, Biosafety in Microbiological and Biomedical Laboratories, 2nd Ed., HHS Publication No. 8808395, 1988; 4-6.
4. American Academy of Pediatrics. Summaries of Infectious Disease. In: Pickering LK, ed. 2000 Red Book: Report of the Committee on Infectious Diseases. 25th ed. Elk Grove Village IL: American Academy of Pediatrics; 2000 528.
5. Lauer BA, Reller LB, and Mirrett S, Effect of atmosphere and duration of incubation on primary isolation of group A Streptococci from throat cultures, J. Clin. Microb., 1983; 17:338-340.
6. Wannamaker LW, Differences between streptococcal infections of the throat and of the skin, N. Eng. J. Med., 1970 282: 23-31, 78-85.

ImmunoCard STAT! is a registered trademark of Meridian Bioscience, Inc.

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Licensed under U.S. Patent Nos. 5,714,389; 5,989,921; 6,485,982 and 6,979,576 and related non-U.S. patents and patent applications.

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INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key Guide to symbols (Guida ai simboli, Guide des symbols, Guia de símbolos, Erläuterung der graphischen symbole)

Manufacturer/Fabbricante/Fabricant/ Fabricante/Hersteller	Date of manufacture/Data di fabbricazione/Date de fabrication/Fecha de fabricación/Herstellungsdatum
In vitro diagnostic medical device/ Dispositivo medico-diagnostico in vitro/ Dispositif médical de diagnostic in vitro/Dispositivo médico para diagnóstico in vitro/ In-Vitro-Diagnostikum	Consult Instructions for Use/Consultare le istruzioni per l'uso/Consulter les instructions d'utilisation/Consulte las instrucciones de uso/ Gebrauchsanweisung beachten
Catalogue number/Numero di catalogo/ Référence du catalogue/Numero de catálogo/Bestellnummer	Use By/Utilizzare entro/Utiliser jusque/Fecha de caducidad/Verwendbar bis
For IVD Performance Evaluation Only/Soltanto per valutazione delle prestazioni/Réactifs IVD réservés à l'évaluation des performances/Sólo para evaluación del funcionamiento/Nur zur IVD Leistungsbewertung	Temperature limitation/Limiti di temperatura/ Limites de température/Limite de temmpérature/ Temperaturbegrenzung
Serial number/Numero di serie/Numéro de série/ Número de serie/Seriennummer	Test Device/Dispositivo test/Dispositif de test/ Dispositivo de Prueba/Testgerät
Batch Code/Codice del lotto/Code du lot/Código de lote/chargenbezeichnung	Buffer/Soluzione tampone/Solution tamponnée/Tampón/Puffer
Conjugate/Coniugato/Conjugué/ Conjugado/Konjugat	Substrate/Substrato/Substrat/ Substrato/Substrat
Contains sufficient for <n> tests/Contenuto sufficiente per "n" saggi/Contenu suffisant pour "n" tests/Contenido suficiente para <n> ensayos/ Inhalt ausreichend für <n> Prüfungen	CAUTION: Risk of Danger/ATTENZIONE: Pericolo/ AVERTISSEMENT: Risques de danger/Precaución: Peligroso/ WARNING: Risikogefahr
Single Use Only/Prodotto Monouso/A usage unique/Para Un Solo Uso/nur für die einmalige Anwendung	Do not freeze/Non congelare/Ne pas congeler/ No congelar/Nicht Eingrieren

TUBE Empty Tube/Provetta vuota/Tube vide/Tubo vacío/Leeres Gefäß	CONTROL Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF WASH 20X Wash Buffer Concentration: 20X/ Soluzione di lavaggio 20X/Solution de lavage concentrée 20X/Solución tampón de lavado 20X/20fach konzentriertes Waschkonzentrat	 Do not use if package is damaged/Non utilizzare se la confezione è danneggiata/ne pas utiliser si le paquet est endommagé/No use si el paquete está dañado/ Nicht verwenden, wenn die Verpackung beschädigt ist
CONTROL - Negative control/Controllo negativo/Contrôle négatif/ Control negativo/ Negative Kontrolle	CONTROL + Positive control/Controllo positivo/Contrôle positif/ Control positivo/Positive Kontrolle
BUF RXN Reaction Buffer/Tampone di reazione/ Solution de réaction tamponnée/ Tampón de Reacción/Reaktionspuffer	Rx Only Prescription Use Only/Per l'uso su prescrizione medica/Uniquement sur prescription/Solo Para Uso Por Receta/ verschreibungspflichtig
SOLN STOP Stopping Solution/Soluzione di Stop/ Solution d'arrêt/Solución de parada/ Stopplösung	CONJ ENZ Enzyme Conjugate/Coniugato enzimatico/ Conjugué enzymatique/ Conjugado enzimático/ Enzymkonjugat
CONTROL Assay Control/Controllo del test/ Test de contrôle/Control de Ensayo/Kontrolltest	REAG Reagent/Reagente/Réactifs/ Reactivos/Reagenzien
BUF WASH Wash Buffer/Soluzione di lavaggio/ Solution de lavage/ Tampón de lavado/ Waschkonzentrat	 Warning/Avvertenze/Mise En Garde/ Advertencia/Warnhinweise
DIL SPE Specimen Diluent (or Sample Diluent)/Diluyente del Campione/ Diluant échantillons/Diluyente de muestra/ Probenverdünnungspuffer	DET REAG Detection Reagent/Reagente Diretto/Réactif de Detection/ Reactivo de Detección/Nachweis Reagenz
CH REP Swiss Authorized Representative/ Mandatario svizzero/Mandataire Suisse/ Representante Autorizado Suizo/Schweizer Bevollmächtigter	 Exclamation mark/Punto esclamativo/Point d'exclamation/ Signo de exclamación/ Ausrufezeichen
EC REP Authorized Representative in the European Community/Rappresentante Autorizzato nella Comunità Europea/Mandataire dans la Communauté européenne/ Representante autorizado en la Comunidad Europea/Bevollmächtigter in der Europäischen Gemeinschaft	
CE This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	
SMP PREP DIL SPE Sample Preparation Apparatus containing Sample Diluent/Dispositivo per la preparazione del campione contenente il diluente del campione/Système pour la préparation de l'échantillon, diluant inclus/Aparato para Preparación de Muestra con Diluyente de Muestra/System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet	

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.

MANUFACTURED BY:

SEKISUI Diagnostics, LLC, 6659 Top Gun Street, San Diego, CA 92121 USA Tel: 781-652-7800

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ImmunoCard STAT! è un marchio registrato di Meridian Bioscience, Inc.
ImmunoCard STAT! est une marque déposée de Meridian Bioscience, Inc.
ImmunoCard STAT! es una marca comercial registrada de Meridian Bioscience, Inc.
ImmunoCard STAT! ist ein eingetragenes Warenzeichen von Meridian Bioscience, Inc.