Para-Pak® Modified (Cu) PVA/10% Formalin Systems

REF 300812

IVD

R_x Only

INTENDED USE

Para-Pak PVA based systems provide standardized procedures for the routine collection, transportation, preservation, and examination of stool specimens for intestinal parasites. Kit systems are designed for easy use by individuals not trained in microbiological procedures and afford an excellent means of minimizing the adverse effects of delay in specimen transportation.

SUMMARY AND EXPLANATION OF THE TEST

Diagnosis of intestinal parasitic disease is confirmed by recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts in the clinical parasitology laboratory. Timely collection and transportation of "fresh" stool specimens to the laboratory cannot always be insured. Workload conditions and priorities in clinical laboratories frequently do not permit immediate examination of "fresh" specimens. Procedures such as incubation, refrigeration, or freezing of stool specimens will not guarantee the recovery of all diagnostic stages of all parasites.^{5, 6, 8, 9, 11, 12}

In 1949, Brooke and Goldman described the use of PVA-fixative for the preservation of intestinal protozoa. In 1981, Horen described substituting cupric sulfate for the more toxic mercuric chloride in the preparation of conventional Schaudin's fixative for use in preparing PVA fixative.¹⁰ The efficacy of the two vial PVA/10% Formalin method in the transportation and preservation of stool specimens has since been confirmed by numerous authors.^{2-6, 8, 9}

Proper use of the Para-Pak Systems thus assures the parasitologist that diagnostic stages of intestinal parasites, if present, will be preserved.

BIOLOGICAL PRINCIPLES

Polyvinyl alcohol is a water soluble plastic which when combined with Schaudin's modified fixative (cupric sulfate) provides an acceptable preservative-fixative for protozoan trophozoites. The resulting permanent slide lends itself to commonly used staining procedures such as Wheatley's trichrome or iron hematoxylin.^{2-4, 13} Ten percent formalin preserved specimens may be examined directly or concentrated for recovery of eggs, larvae, and protozoan cysts.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. PVA fixative (item# 9005)
2. 10% Formalin (item# 9004)

Each kit consists of one vial containing Modified PVA fixative and one vial containing 10% Formalin preservative. Single vial cases are also available. Simple directions for patients and nursing personnel are also provided.

MATERIALS NOT PROVIDED

1. Ethyl acetate (suggested) or ether (optional)
2. Zinc sulfate solution (specific gravity = 1.18)
3. Physiological saline
4. Cotton tipped applicator sticks
5. Microscope slides and coverslips
6. Centrifuge
7. Microscope
8. Transfer pipettes

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Avoid contact of fixative solutions with the skin and eyes. Should contact occur, flush with running water. If irritation should develop, see a physician.
3. Fixative solutions are poisonous. If ingested, dilute by drinking milk or water. Then call local poison center or physician immediately.
4. Due to the infectious nature of unpreserved stools, care and handwashing should be employed when the specimen is collected and manipulated.
5. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800-343-3858.
6. IMPORTANT: See SDS for additional safety and hazard information.

HAZARD and PRECAUTIONARY STATEMENTS

 <p>Para-Pak 10% Formalin</p>	<p>Signal Word Danger</p> <p>Hazard Statements H301 - Toxic if swallowed H311 - Toxic in contact with skin H317 - May cause an allergic skin reaction H330 - Fatal if inhaled H331 - Toxic if inhaled H341 - Suspected of causing genetic defects H350 - May cause cancer H315 - Causes skin irritation H318 - Causes serious eye damage H370 - Causes damage to organs H402 - Harmful to aquatic life Contains Formaldehyde, Methyl alcohol</p> <p>Precautionary Statements - EU (§28, 1272/2008) P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor. P321 - Specific treatment (see Section 4 and Section 11 on this label). P280 - Wear protective gloves/protective clothing/eye protection/face protection. P321 - Specific treatment (see supplemental instructions on the administration of antidotes on this label). P361 + P364 - Take off immediately all contaminated clothing and wash it before reuse. P260 - Do not breathe dust/fume/gas/mist/vapors/spray. P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing. P310 - Immediately call a POISON CENTER or doctor. P320 - Specific treatment is urgent (see Section 4 and Section 11 on this label). P403 + P233 - Store in a well-ventilated place. Keep container tightly closed. P280 - Wear eye protection/ face protection. P201 - Obtain special instructions before use. P308 + P313 - IF exposed or concerned: Get medical advice/attention. P281 - Use personal protective equipment as required. P202 - Do not handle until all safety precautions have been read and understood. P270 - Do not eat, drink or smoke when using this product. P271 - Use only outdoors or in a well-ventilated area. P284 - Wear respiratory protection. P272 - Contaminated work clothing should not be allowed out of the workplace. P307 + P311 - IF exposed: Call a POISON CENTER or doctor/physician. P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352 - IF ON SKIN: Wash with plenty of soap and water. P312 - Call a POISON CENTER or doctor/physician if you feel unwell. P361 - Remove/Take off immediately all contaminated clothing. P332 + P313 - If skin irritation occurs: Get medical advice/attention. P363 - Wash contaminated clothing before reuse. P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. P405 - Store locked up. P501 - Dispose of contents/ container to an approved waste disposal plant.</p>
--	---

 <p>Para-Pak Modified PVA (Cu) Fixative</p>	<p>Signal word Danger</p> <p>Hazard statements H225 - Highly flammable liquid and vapor H319 - Causes serious eye irritation H336 - May cause drowsiness or dizziness H411 - Toxic to aquatic life with long lasting effects Contains Isopropyl alcohol</p> <p>Precautionary Statements - EU (§28, 1272/2008) P264 - Wash face, hands and any exposed skin thoroughly after handling. P280 - Wear eye protection/ face protection. P337 + P313 - If eye irritation persists: Get medical advice/attention. P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P370 + P378 - In case of fire: Use CO2, dry chemical, or foam to extinguish.</p>
--	--

SHELF LIFE AND STORAGE

Shelf life of the Para-Pak System is indicated on the outer package label. Store at room temperature (15-30 C). Excessive heat or cold should be avoided.

SPECIMEN COLLECTION AND PREPARATION

1. The patient should be cautioned against the use of antacids, barium, bismuth, antidiarrheal medication, or oily laxatives prior to collection of the specimen.
2. To assure recovery of parasitic elements which are passed intermittently and in fluctuating numbers, three specimens spaced a few days apart must be examined. In the case of hospitalized patients it is suggested that all fecal passages be collected for a designated length of time to avoid prolonging the hospital stay.^{5,9}
3. The specimen is ideally passed into a bedpan but must not be contaminated with urine. Alternatively, plastic wrap may be placed in the toilet seat opening and the specimen passed into the bag. A thoroughly cleaned and dried milk carton, cut so as to remove the upper two thirds of the carton may also be used. It will be easier to collect the specimen if the water supply to the toilet is shut off and water drained from the bowl.
4. An appropriate (i.e. bloody, slimy, watery) area of stool should be selected and sampled with the collection spoons provided in the caps of the vials. Sufficient stool is added to each vial to bring the liquid level up to the "Fill to Here" line. This will result in approximately 5 mL of sample. To insure ideal sampling of a formed stool, material should be removed from the sides, ends, and middle of the bolus.
5. Agitate each specimen with the spoon along the sides of the vial, tighten the cap and shake firmly to insure that the specimen is adequately mixed. When mixing is completed the specimen should appear homogeneous.
6. Complete the label on each vial and seal the vials in the plastic bag.

TEST PROCEDURE

The Para-Pak System lends itself to a wide variety of procedures in common use. The following discussion is not exhaustive and alternatives may be found in the literature cited. While variations exist from lab to lab, a thorough examination should include at LEAST four steps:

1. Gross examination: record the presence of blood, worms, mucus, or proglottids.
2. Direct microscopic examination from the 10% Formalin preserved specimen:
 - a. Place a clean glass slide on a sheet of newsprint.
 - b. Add a drop of saline (iodine may be substituted if desired) to the slide.
 - c. Add a representative sample of formalin preserved specimen to the drop of saline and mix thoroughly with the collecting spoon. The newsprint must be just legible through the slide.
 - d. Place a double width coverslip on the suspension and examine immediately.
3. Permanent slides for staining with Wheatley's Trichrome Stain (Meridian Cat. #400101), iron-hematoxylin, etc.:
 - a. Pour some of the Modified PVA-fixed material onto a paper towel and allow to stand for three minutes. This will absorb out excess PVA and is considered critical to obtain the best possible staining.⁸
 - b. Using an applicator stick or brush, spread (avoid smearing) some of the specimen from the paper towel onto one or more clean glass slides. For best adherence, spread the material to the edge of the slide.
 - c. The slides are dried overnight at room temperature or for several hours in a 37 C incubator or slide warmer. Accelerated drying may cause some morphological distortion. The slides must be dried thoroughly to avoid washing the film off during staining.
 - d. Once slides prepared in Modified PVA have dried completely, staining may begin by placing the slides directly into the trichrome stain and proceeding as normal.

NOTE: Slides made from very watery specimens may require an additional day to dry completely.
4. Concentration procedures: One or more concentration procedures should be employed. No one concentration procedure works equally well for all parasites,^{1, 5, 11} however, two in common usage that lend themselves well with the Para-Pak System are:
 - A. Formalin-ether (ethyl acetate^{7, 14}) sedimentation:
 1. Mix the 10% Formalin or Modified PVA specimen thoroughly. The specimen is now ready for processing with the Para-Pak® Macro-Con® or CON-Trate® Stool Concentration Systems. See the appropriate package insert for further directions. If Para-Pak Macro-Con or CON-Trate is not available, a sufficient quantity of specimen must be strained into a 15 mL conical centrifuge tube through one layer of narrow mesh or two layers of wide mesh gauze to provide the amount of sediment required in Step No. 2. This amount will vary with the size and density of the specimen.
 2. Add saline, mix thoroughly and centrifuge at 500 xg (1800-200 rpm for most tabletop centrifuges) for 10 minutes. If the resulting sediment is not approximately 1 mL, resuspend, add or remove specimen and recentrifuge.
 3. Decant the supernatant fluid. A second wash may be used if desired.
 4. Add approximately 10 mL of 10% Formalin, mix thoroughly and allow to stand five minutes.
 5. Add 3 mL of ethyl acetate or ether then stopper and shake vigorously for at least 30 seconds. Carefully remove the stopper.
 6. Centrifuge at 500 xg for 10 minutes (1800-2200 rpm).
 7. Four layers will be apparent:
 - a. Top layer: ethyl acetate or ether
 - b. Second layer: plug of debris
 - c. Third layer: formalin
 - d. Bottom layer: sediment
 8. After ringing the plug of debris from the sides of the tube with an applicator stick, carefully decant the top three layers. While keeping the tube inverted, a cotton swab may be used to remove debris sticking to the sides of the tube. This is particularly important for obtaining suitable results with ethyl acetate and avoids solvent bubbles in the wet mount.
 9. Add a few drops of physiological saline or 10% Formalin to resuspend the remaining sediment. If the resulting slides are too dense (newsprint should be legible through them) more saline or formalin may be added.
 10. Iodine and saline mounts are suggested for microscopic examination.
 - B. Zinc Sulfate flotation:
 1. Thoroughly mix a representative portion of the 10% Formalin stool suspension or fresh unpreserved specimen in a 15 mL centrifuge tube and q.s. with tap water to approximately 10-12 mL. The amount of specimen to use will vary with its size and density.
 2. Centrifuge one minute at 1000-1200 xg.
 3. If sediment is about one milliliter in volume, decant the supernatant fluid. Otherwise adjust the density of the suspension by adding material from the 10% Formalin suspension or diluting with more water. If adjustment of the sediment was necessary, or if the stool is very oily, repeat the wash procedure.
 4. When using formalinized specimens, the specific gravity of the zinc sulfate solution must be adjusted to 1.2.^{1, 5} Fill the tube about half full with zinc sulfate solution and resuspend the sediment by mixing thoroughly with applicator sticks.
 5. Add additional zinc sulfate solution to within one inch of the top.
 6. Centrifuge one minute at 1000-1200 xg.
 7. Carefully remove the tube from the centrifuge and, avoiding agitation, place it upright in a test tube rack or other suitable holder.
 8. Carefully fill the tube to the brim with zinc sulfate solution. Do not allow any overflow.
 9. A clean coverslip may now be placed on top of the tube. If the coverslip does not contact the meniscus of the liquid, carefully add more zinc sulfate solution until it does.
 10. Do not disturb the tube or coverslip for ten minutes.
 11. With a quick, deft motion remove the coverslip straight upward so that a drop of liquid containing eggs and cysts adheres to the center of the coverslip.
 12. The coverslip may now be placed on a clean glass slide. If an iodine mount is desired, place a small drop of iodine on the slide prior to adding the coverslip. Sealing the edge of the slide with Vaspar (vaseline/paraffin (1:1) mixture) will prevent drying and distortion of the larger eggs.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Visual inspection: Vials of Modified PVA should contain approximately 15 mL of clear blue fluid, to insure 1:3 stool to preservative ratio.
2. If gelled, the fixative may be liquified by placing in a 50 C water bath until clear and fluid.
3. On prolonged storage or exposure to cold temperatures, a small white or blue precipitate may form. This will not affect adhesive, staining, or fixative properties. The precipitate may go back into solution if placed in a 42 C incubator for 2-3 days without affecting the performance of the product.
4. When a PVA fixed film of stock trophozoite or human buffy coat is stained, the organisms or cells should appear well fixed and defined.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

Para-Pak® Modified (Cu) PVA/10% Formalin Systems

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REFERENCES

1. Bartlett, Marilyn, et al. Comparative evaluation of a modified zinc sulfate flotation technique. *J. Clin Microbiol* 1978;7:524-528.
2. Brook, MM. and M. Goldman. Polyvinyl alcohol fixative as a preservative and adhesive for protozoa in dysenteric stools, and other liquid materials. *J Lab and Clin Med* 1944;34:1554-1560.
3. Brooke, MM. and C. Norman. The effectiveness of the PVA fixative technique in revealing intestinal amebae in diagnostic cultures. *Am J Trop Med Hyg* 1955;4:479-482.
4. Brooke, MM. PVA fixative technique in laboratory confirmation of amebiasis. *Triangle* 1960;4:326-335.
5. Brooke, MM. "Intestinal and Urogenital Protozoa" in *Manual of Clinical Microbiology*, ASM, Washington, D.C., Second Edition, 1974:582-601.
6. Burrows, RB. *Microscopic Diagnosis of the Parasites of Man*, Yale University Press, New Haven. 1965.
7. Erdman, Dean. Clinical comparison of ethyl acetate and diethyl ether in the formalin-ether sedimentation technique. *J. Clin Microbiol* 1981;14:483-485.
8. Garcia, LS. and L. Ash. *Diagnostic Parasitology*, C. V. Mosby Co., St. Louis, 1979:pp 9-24.
9. Harper, K., et al. Advantages of the PVA fixative two-bottle stool collection technique in the detection and identification of intestinal parasites. *Pub Health Lab* 1957;15:96-108.
10. Horen, WP. Modification of Schaudin's fixative. *J Clin Microbiol* 1981;13:204-205.
11. Scholten, T. H. and J. Yang. Evaluation of unpreserved and preserved stools for the detection and identification of intestinal parasites. *Am J Clin Path* 1974;62:563-567.
12. Swartzwalder, J. Clyde, GW. Hunter, WW. Frye. *Manual of Tropical Medicine*, W. B. Saunders Co., Philadelphia. 1966.
13. Simitch, T. and Z. Petrovich. Longevite de la forme vegetive de dysenterioe dans divers mileaux. *Arch Inst Pasteur d'Algerie* 1953;31:375-380.
14. Young, Kirk H., et al. Ethyl acetate as a substitute for diethyl ether in the formalin-ether sedimentation technique. *J Clin Microbiol*. 1979;10:852-853.

SN10710

REV. 08/25

Manufactured By

Meridian Bioscience, Inc.

3471 River Hills Drive
Cincinnati, OHIO - 45244 USA
www.meridianbioscience.com


Contacts:

Main Telephone (+1) 513.271.3700
Customer Service/Orders 800.543.1980
Technical Support Center 800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124
E-mail: info@meridianbioscience.com

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/package of this product:

Key guide to symbols (Guida ai simboli/Guide des symboles, Guia de simbolos, Erläuterung der graphischen symbole)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.