Meridian Viral Transport

REF 505021

IVD

Rx Only

INTENDED USE

Meridian Viral Transport is a system which contains all materials necessary to collect and transport clinical specimens for viral isolation.

SUMMARY AND EXPLANATION OF THE TEST

Classically, the method of choice for the demonstration of the presence of virus in patient specimens involves tissue culture isolation. The efficacy of this method depends upon the collection of a specimen containing infectious virus and the maintenance of the infectivity during transport to the testing laboratory.

Many viruses (such as respiratory syncytial virus⁵) are labile and specimens for isolation require a transport media to stabilize their infectivity.¹⁻⁶

BIOLOGICAL PRINCIPLES

Meridian Viral Transport contains a sterile swab for collection of a clinical specimen and a tube containing a stabilizing medium which has been proven effective in sustaining viral infectivity.¹

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

Meridian Viral Transport contains all materials needed for collection of routine swab-type specimens.

Contents

1. Sterile swab – Flocked polyester head on a plastic shaft
2. Transport Tube – containing 1.5 mL of medium
Medium formulation per liter:

Tryptose phosphate broth	29.5 g/L
Gelatin	5.0 g/L
Gentamycin sulfate	50.0 mg/L

MATERIALS NOT PROVIDED

Certain types of specimens (eg, lesion scrapings, vesicular fluid) may require the use of additional materials (scalpel, hypodermic needle). The materials that will be required depend on the type of the specimen and the methodology of the physician.

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Meridian Viral Transport tubes containing patient specimens should be considered potentially infectious and handled accordingly. All clinical specimens should be autoclaved before disposal.

HAZARD and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

Store at 2-25 C before use. Inspect each tube visually before use. DO NOT USE tubes showing evidence of contamination (cloudiness or precipitate). Do not use beyond the labeled expiration date on the outer box label.

TEST PROCEDURE

Specimen Collection

1. Peel open the package and remove the sterile swab. Collect a patient specimen with the swab employing standard clinical procedures.¹⁻⁶
2. Remove the cap from the tube. Insert the inoculated swab and break off the swab shaft at score mark by bending shaft against the tube rim.
3. Replace the cap on the tube and seat it securely.
4. Label tube with patient identification and forward it to the laboratory for testing.

Specimen Transport

The mode of transport of patient specimens to the testing laboratory depends upon the virus being sought. The following is intended as a general guideline applicable to Herpes Simplex, as well as a number of other clinically relevant viruses. (See also section on LIMITATIONS)

1. Specimens should be placed promptly in Meridian Viral Transport for best results.
2. Specimens should be held at refrigerator temperature (2-8 C) until they reach the testing laboratory.
3. If specimens are to be delayed longer than 72 hours, they should be frozen and shipped at ≤ -70 C. (See section on LIMITATIONS for further information regarding storage.)
Note: Some virus infectivity will be lost upon freezing. Avoid repeated freezing and thawing or storage in self-defrosting freezers.
4. Leftover specimens can be stored frozen at ≤ -70 C for future reference.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

The preservation capabilities of each lot of Meridian Viral Transport are confirmed before release in tests for Herpes Simplex Virus (HSV) infectivity. For most applications, no additional controls need be performed by the user.

LIMITATIONS OF THE PROCEDURE

1. The transport and storage of many viruses requires special handling. The user should consult the standard references listed in the References for additional information regarding proper transportation conditions.
2. Viral transport media may be employed for collecting specimens for viral detection assays. Viral detection assays may have more stringent collection, storage and handling requirements than defined in this insert. In such cases, follow the more stringent criteria of the immunoassay. This product was cleared for use with TRU FLU and TRU RSV rapid immunoassays. Consult the package insert for those products for optimum sample use, storage and transport conditions.

EXPECTED VALUES AND SPECIFIC PERFORMANCE CHARACTERISTICS

The tryptose phosphate broth contained in Meridian Viral Transport tubes is a standard broth-based medium for the transport of clinical specimens intended for virus isolation procedures. (See Reference 1.)

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REFERENCIAS

- Johnson FB. Transport of viral specimens. Clin Microbiol Rev 1990; 3:120-131.
- Smith TF. Specimen requirements: Selection, collection, transport and processing. In: Specter S, Hodinka RL, Young SA, Eds. Clinical virology manual, ed 3, Washington DC, ASM Press: 2000; pp11-26.
- Forman MS, Valsamakis A. Specimen collection, transport and processing: virology. In: Murray PR et al. Eds. Manual of Clinical microbiology, ed 9. Washington DC: ASM Press, 2007; pp 1284-1296.
- Atmar RL. Influenza viruses. In: Murray PR et al. Eds. Manual of Clinical microbiology, ed 9. Washington DC: ASM Press, 2007; pp 1340-1351.
- Tang Y-M, Crowe JE Jr. Respiratory syncytial virus and human metapneumovirus. In: Murray PR et al. Eds. Manual of Clinical microbiology, ed 9. Washington DC: ASM Press, 2007; pp 1361-1377.
- Jerome KR, Morrow RA. Herpes simplex viruses and herpes B virus. In: Murray PR et al. Eds. Manual of Clinical microbiology, ed 9. Washington DC: ASM Press, 2007; pp 1523-1531.

SN10867

REV. 09/22

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guía de símbolos, Guide des symboles, Guia de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call customer Service Department at 800-543-1980.

 Manufactured By	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
----------------------------	---