immunocard STAT![®]

E. coli O157 Plus

A Rapid Immunoassay for the Detection of Shiga Toxin-Producing *E. coli* O157 in Human Stool Specimens and Culture

REF 750530

IVD

Rx Only

INTENDED USE
The Immunocard STAT! *E. coli* O157 Plus is a rapid test for the detection of antigens from Shiga toxin-producing *E. coli* O157 as an aid in the diagnosis of *E. coli* O157:H7 infection. The test can be used to directly test stool specimens, stool in modified Cary-Blair medium or confirmatory stool cultures grown in MacConkey broth of sorbitol MacConkey (SMAC) plates.

In the US, Immunocard STAT! *E. coli* O157 Plus is not intended for point-of-care use. This device is intended for moderately complex laboratories. In Canada, this device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

Escherichia coli is a common microorganism found in normal human flora. While many different serotypes have been identified, certain strains have acquired genes for toxins from other bacterial species via bacteriophage plasmids. Shiga toxin-producing *E. coli* O157 usually carries one or more genes for Shiga toxin, placing it in a category of pathogens known as Shiga toxin-producing enterohemorrhagic *E. coli* (STEC). STEC strains of *E. coli* cause a wide spectrum of human diseases, including bloody and non-bloody diarrhea, hemorrhagic colitis, kidney failure, hemolytic uremic syndrome (HUS), and death.¹⁻⁵ Morbidity and mortality rates in children, infants and the elderly are especially high. In addition, many states and institutions require a negative test result prior to readmitting infected children to schools or day-care facilities. Outbreaks attributed to *E. coli* O157 have been associated most commonly with undercooked beef and dairy products, but have also been traced to leaf lettuce, radish sprouts, alfalfa sprouts, unpasteurized apple cider⁶⁻⁸ and contaminated water.^{9, 10}

To date, *E. coli* O157:H7 has been the most commonly described STEC serotype.^{6, 11, 12} The Centers for Disease Control (CDC) recommends that all bloody stools submitted for routine stool culture be examined for the presence of *E. coli* O157:H7.⁶ Some investigators recommend that all stools submitted for routine stool culture be examined for *E. coli* O157:H7.⁴ Most sorbitol-nonfermenting *E. coli* O157 are motile and possess the H7 antigen. However, motility can be difficult to elicit, resulting in failure to detect the H7 antigen. Occasionally, a sorbitol-nonfermenting strain of *E. coli* O157 is either nonmotile (O157:NM), or motile with a nontypable flagellar antigen (O157:H-). Clinical laboratory evaluation of a strain need not include H7 antigen determination.⁴ Sorbitol-nonfermenting *E. coli* O157 should be considered pathogenic, whether or not the H7 antigen is detected.⁴ Determination of H7 may be helpful epidemiologically, or for state and federal reporting purposes.

Current methods for detection of Shiga toxin-producing *E. coli* O157 include 1) culture on SMAC plates followed by confirmatory testing, 2) enzyme immunoassays for either O157 antigens or toxins, and 3) broth amplification followed by either EIA, latex, or subculture. The EIA methods require approximately one hour to complete and require a high skill level with multiple process steps. The culture methods required at least one overnight incubation and skill in identifying sorbitol negative colonies as well as biochemical and/or immunological identification. However, culture methods amplify the number of organisms present to detectable levels in specimens with low numbers of *E. coli* O157. Immunocard STAT! *E. coli* O157 Plus was designed for the rapid and accurate detection of Shiga toxin-producing *E. coli* O157. Accurate and rapid detection of Shiga toxin-producing *E. coli* O157 is necessary to prevent further spread of the organism, avoid unnecessary procedures such as surgery and to guide treatment decisions.

BIOLOGICAL PRINCIPLES

Stool or culture material are prepared / diluted and added to the sample port of the device. The sample mobilizes gold particles, coated with monoclonal antibody specific for the *E. coli* O157 lipopolysaccharide, and migrates along the membrane through the Test and Control zones. The Test zone contains immobilized monoclonal antibody specific for an epitope common to Shiga toxin-producing *E. coli*. After 10 minutes the Test and Control zones are observed for the presence of red/purple lines across the membrane surface. If a Shiga toxin-producing *E. coli* O157 is present in the sample, a complex is formed between the capture antibody, the Shiga toxin-producing *E. coli* O157, and the monoclonal antibody-gold conjugate which can be seen visually as a red/purple line in the Test zone. No red/purple line in the Test zone indicates a negative result. The Control line serves as a procedural control, to assure that the sample has migrated the appropriate distance along the membrane.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Immunocard STAT! *E. coli* O157 Plus Devices** – Individually foil pouched devices containing immobilized monoclonal antibody specific for an epitope common to Shiga toxin-producing *E. coli* (Test zone), goat anti-mouse IgG antibody (Control zone) and monoclonal antibody, specific for the *E. coli* O157 lipopolysaccharide, conjugated to gold particles.
2. **Positive Control** – Inactivated *E. coli* O157:H7 in a solution containing 0.094% sodium azide as a preservative.
3. **Negative Control** – Inactivated *E. coli* O157:H12 in a solution containing 0.094% sodium azide as a preservative.
4. **Sample Diluent** – Solution containing 0.094% sodium azide as a preservative.
5. **Transfer pipettes**

MATERIALS NOT PROVIDED

1. 12 x 75 mm test tubes
2. Timer

OPTIONAL MATERIALS NOT PROVIDED

1. MacConkey Broth (5.0 mL/tube)
2. MacConkey Agar w/Sorbitol (SMAC)
3. Incubator
4. Swabs
5. McFarland Standards (#2 through #4)
6. Modified Cary-Blair medium

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Directions should be read and followed carefully.
3. Reagent concentration, incubation times and temperatures (21–27 C) have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications. Once the assay has been started, complete all subsequent steps without interruption.
4. All reagents should be gently mixed and at 21–27 C before use.
5. Hold Positive and Negative Control vials vertically at suitable distance above the well to insure proper drop size and delivery.
6. Do not allow the tips of the vial or pipette to touch the Sample Port.
7. All reagents are provided already diluted to the proper concentration. Do not dilute further.
8. Pipettes that are supplied with the kit should be marked at the intervals shown in the pipette diagram in PROCEDURAL NOTES. When instructed to use the pipette supplied with the kit, do not use transfer pipettes with markings that differ from the diagram. Use one per specimen.
9. Do not interchange reagents from different lot numbers.
10. Test devices are packaged in foil pouches that exclude moisture during storage. Inspect each foil pouch before opening. Do not use test devices from pouches that have holes in the foil or where the pouch has not been completely sealed. False negative reactions may result if test components and reagents are improperly stored.
11. Positive and Negative Control reagents contain inactivated *E. coli*. However, they should be handled as potential biohazards.
12. Do not use kit or components beyond their assigned expiration dates.
13. Do not use vials that lack a label, a lot number, or an expiration date.
14. Replace colored caps on correct vials.
15. Patient specimens and used Immunocard STAT! *E. coli* O157 Plus devices may contain infectious agents and should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories."
16. Some reagents contain the preservative sodium azide which is a skin irritant. Avoid contact with reagents. Disposal of reagents containing sodium azide into drains consisting of lead or copper plumbing can result in the formation of explosive metal oxides. Eliminate the build-up of oxides by flushing drains with large volumes of water during disposal.
17. Dispose of stool specimens as potentially biohazardous materials. Inactivate cultures by autoclaving for a minimum of 15 minutes at 121 C before disposal.
18. Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
19. Stool must be mixed thoroughly, regardless of consistency, to ensure a representative sample prior to pipetting.

HAZARDS and PRECAUTIONARY STATEMENTS


There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C when not in use.

PROCEDURAL NOTES

- 1. The Immunocard STAT! *E. coli* O157 Plus format is diagrammed below.
- 2. Batch processing of samples or controls is possible provided correct incubation time is maintained for each device.
- 3. The Control zone of each device is a procedural control to assure that the sample has migrated sufficiently in the device to permit a valid test result to be read.

REAGENT PREPARATION

- 1. All reagents come ready to use (no further dilution is required).
- 2. Allow kit components to reach 21-27 C prior to use.
- 3. Gently mix liquid reagents prior to use.

SPECIMEN COLLECTION AND PREPARATION

Stool can be tested either directly, as obtained from modified Cary-Blair medium or after culture in MacConkey broth or on a SMAC plate. Stool samples should be received in an airtight container and may be stored up to 72 hours at 2-8 C prior to testing or culture. Freeze specimens immediately upon receipt if testing cannot be completed within 72 hours. For culture confirmation, MacConkey broth or a SMAC plate should be set up as soon as possible after receiving a stool specimen. Freezing specimens prior to culture may affect recovery of viable organism.


- 1. **BROTH ENRICHMENT:** Add approximately 50 µL of stool to 5 mL MacConkey Broth and incubate at 42 C ± 1 C for 16-24 hours.
- 2. **PLATE CULTURE:** Streak stool, broth or culture isolate onto a SMAC plate and incubate at 35 C ± 1 C for 16-24 hours.

TEST PROCEDURE

- A. DIRECT STOOL:**
- 1. Add 350 µL of Sample Diluent to a test tube using the dropper assembly.
 - 2. Mix stool as thoroughly as possible prior to pipetting:
 - a. For liquid or semi-solid stools draw the mixed stool specimen to the first mark (25 µL) on a transfer pipette. Add to Sample Diluent and mix. Note: Do not pipet more than 25 µL of stool. Over-inoculation with stool may produce invalid results. Leave the transfer pipette in the tube for future use.
 - b. For non-pipetable stools, use a wooden applicator stick and transfer a small portion (2 mm diameter) of thoroughly mixed stool into the Sample Diluent and emulsify the stool thoroughly with the applicator stick. Place a transfer pipette in the tube for future use.
 - 3. Open a foil pouch and remove the device. Vortex diluted specimen for 10 seconds, then add 150 µL (second mark on transfer pipette) of diluted stool to the sample port of the device.
 - 4. Incubate for 10 minutes at 21-27 C. Note: during the 10 minute incubation, diluted specimen must move past the Control line.
 - 5. Visually read the Control and Test zones for the presence or absence of a red/purple line at the end of the incubation period.
- B. BROTH ENRICHMENT or CARY-BLAIR SPECIMEN:**
- 1. Add 350 µL Sample Diluent to a test tube using the dropper assembly.
 - 2. Draw broth culture to the base of the bulb of a transfer pipette (approximately 350 µL) and add it to the Sample Diluent and mix. For Cary-Blair specimens, draw inoculated Cary-Blair medium to the second line of the transfer pipette (approximately 150 µL) and add it to the Sample Diluent and mix. Leave the transfer pipette in the tube for future use.
 - 3. Open a foil pouch and remove the device. Vortex the diluted specimen for 10 seconds, then add 150 µL (second mark on transfer pipette) of the diluted specimen to the sample port of the device.
 - 4. Incubate for 10 minutes at 21-27 C. Note: During the 10 minute incubation, the diluted specimen must move past the Control line.
 - 5. Visually read the Control and Test zones for the presence or absence of a red/purple line at the end of the incubation period.
- C. PLATE CULTURE:**
- 1. Add 700 µL (2 x 350 µL) of Sample Diluent to a test tube using the dropper assembly.
 - 2. Use a swab to make a heavy suspension (McFarland Standard 2-4) of suspected growth from the SMAC plate in the Sample Diluent. Mix gently.
 - 3. Open a foil pouch and remove the device. Mix the diluted specimen gently, then add 150 µL (second mark on transfer pipette) to sample port of device.
 - 4. Incubate for 10 minutes at 21-27 C. Note: During the 10 minute incubation, the diluted specimen must move past the Control line.
 - 5. Visually read the Control and Test zones for the presence or absence of a red/purple line at the end of the incubation period.

INTERPRETATION OF RESULTS

NOTE: Test results should be read immediately following the 10 minutes incubation period. Results may not be read after an additional 10 minutes.

Positive Test Result: Visually detectable red/purple Test and Control lines. A positive result indicates the presence of antigens from Shiga toxin-producing *E. coli* O157.

NOTE: In view of the epidemiological importance of obtaining bacterial isolates it is recommended that all toxin positive samples be submitted for isolation of Shiga toxin positive organisms. We suggest that individual microbiology laboratories coordinate bacterial isolation with their local state health laboratories.

Negative Test Result: Visually detectable red/purple Control line. No red/purple Test line present. A negative result indicates that antigens from Shiga toxin-producing *E. coli* O157 are absent or below the level of detection

Invalid Test Result: No visually detectable red/purple Control line, with or without a visually detectable red/purple Test line. Invalid test results may be due to a Reagent/Device problem, a procedural error (repeat the test), or over-inoculation of stool into Sample Diluent (re-dilute stool and repeat the test).

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

The Positive and Negative Controls should be assayed in accordance with local, state, and federal regulations. Your laboratory may have established control testing requirements beyond this recommendation.

Although specimen matrix interference has not been observed in this assay, the control reagents furnished are in an aqueous solution matrix, and may not adequately control for the specimen matrix. For any diagnostic assay, the National Committee for Clinical Laboratory Standards, guideline C24-A, indicates that matrix control materials must be used when available. To comply with this guideline, it may be necessary for the user to provide control material in specimen matrix.

CONTROL REAGENT TESTING

Add three drops Positive Control or Negative Control directly to bottom port of separate devices.

- 1. The Positive Control should yield visually detectable red/purple Test and Control lines after 10 minutes incubation at 21-27 C.
 - 2. The Negative Control should yield a visually detectable red/purple Control line and no Test line should be present after 10 minutes incubation at 21-27 C.
- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage.

If the expected reactions are not observed and the reagents are still within their expiration date, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

Shiga toxin-producing *E. coli* O157 causes both sporadic and epidemic outbreaks of disease. While the incidence of disease peaks in the summer months, cases occur year round. Due to the focal nature of the disease, incidence values range widely with geographic locations, from zero to 3-4% in severe outbreaks. For this reason, predictive values for the Immunocard STAT! *E. coli* O157 Plus test are subject to variation. The tables below give the expected predictive values based upon incidence.

O157:H7 Incidence	Stool Specimens		MacConkey Broth		SMAC Plate	
	PPV	NPV	PPV	NPV	PPV	NPV
4.0%	75.5%	99.2%	100.0%	100.0%	100.0%	100.0%
3.0%	69.6%	99.4%	100.0%	100.0%	100.0%	100.0%
2.0%	60.2%	99.6%	100.0%	100.0%	100.0%	100.0%
1.0%	42.8%	99.8%	100.0%	100.0%	100.0%	100.0%
0.5%	27.1%	99.9%	100.0%	100.0%	100.0%	100.0%
0.2%	12.9%	100.0%	100.0%	100.0%	100.0%	100.0%
0.1%	6.9%	100.0%	100.0%	100.0%	100.0%	100.0%

LIMITATIONS OF THE PROCEDURE

- 1. Test results should be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
- 2. The limit of detection was determined utilizing *E. coli* O157:H7 diluted in saline. The limit of detection for *E. coli* O157:H7 ATCC strains 43888 and 43895 was 8.3 x 10⁵ cfu/mL and 7.67 x 10⁵ cfu/mL respectively.
- 3. A negative test result indicates the absence of Shiga toxin-producing *E. coli* O157 antigens, or that the levels of antigens are below that which can be detected by this test. Further testing should be performed with negative specimens to determine the cause of diarrhea. Optimal sensitivity is obtained by culture confirmatory methods.
- 4. Many states require a negative STEC or *E. coli* O157 test result before a child that has experienced a STEC *E. coli* O157 infection may return to school or day care. For this purpose, to reduce the risk of sending infectious children back into these environments, a culture confirmatory method should be utilized.
- 5. Laboratories should comply with all local state and federal regulations for reporting STEC of *E. coli* O157 disease.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Immunocard STAT! *E. coli* O157 Plus test was evaluated at six locations.^{9, 10, 13-16} These sites included four Children's Hospitals and two general hospitals. Immunocard STAT! *E. coli* O157 Plus was tested on direct stool (n=299 fresh, 483 frozen), 24 hour broth cultures (n=277), and on sorbitol negative colonies from SMAC plates (n=52).

In the table below, Immunocard STAT! *E. coli* O157 Plus results are compared to latex testing of sorbitol negative colonies from a SMAC plate.

O157:H7 Culture Results	Immunocard STAT! <i>E. coli</i> O157 Plus Results					
	Direct Stool		MacConkey Broth		SMAC Plate	
	Pos.	Neg.	Pos.	Neg.	Pos.	Neg.
Positive	341	77	17	0	20	0
Negative	5	359	0	261	0	32
Sensitivity	82%	78%-85%	100%	81%-100%	100%	83%-100%
Specificity	99%	97%-100%	100%	99%-100%	100%	89%-100%

The Immunocard STAT! *E. coli* O157 Plus test was 82% sensitive and 99% specific when tested directly on stool specimens. The sensitivity of the assay was greater when culture confirmation methods were used. The sensitivity and specificity of the Immunocard STAT! *E. coli* O157 Plus test were 100% when MacConkey broth and colonies from SMAC plates were tested.

MODIFIED CARY-BLAIR SPECIMEN TESTING

Specimens in modified Cary-Blair transport medium were evaluated with the Immunocard STAT! *E. coli* O157 Plus test. Specimens should be kept at 4 C until tested. Specimens should be tested preferably within 24 hours, but no later than 72 hours.

Immunocard STAT! <i>E. coli</i> O157 Plus	MODIFIED CARY-BLAIR SPECIMEN TESTING CULTURE		
	Positive	Negative	Total
Positive	33	3	36
Negative	6	400	406
Total	39	403	442

SENSITIVITY & SPECIFICITY CARY-BLAIR EVALUATION		
		95% CONFIDENCE INTERVAL CALCULATED USING THE EXACT METHOD
Sensitivity	70%-93%	70%-93%
Specificity	98%-100%	98%-100%

REPRODUCIBILITY

Three physician office laboratories tested six specimens and the two controls, in triplicate, on each of three different days. The Test line of the Immunocard STAT! *E. coli* O157 Plus device showed 100% reproducibility with two negative, two low positive and two medium positive stools. Control line reproducibility was 100% with all specimens.

ASSAY SPECIFICITY

The specificity of Immunocard STAT! *E. coli* O157 Plus was tested by utilizing the following bacterial or viral strains. Positive and negative stools were spiked with ≥ 1 x 10⁸ organisms/mL and tested by Immunocard STAT! *E. coli* O157 Plus. *E. coli* O157:H7 gave a positive result when tested. All other organisms were found to be negative when spiked into the negative stool. In addition, they did not interfere with the positive specimen.¹⁷

Microorganism or virus (# strains tested)

<i>Adenovirus</i> 40 & 41 (2)	<i>Pseudomonas aeruginosa</i> (1)
<i>Campylobacter coli</i> (1)	<i>Pseudomonas fluorescens</i> (2)
<i>Campylobacter fetus</i> (1)	<i>Rotavirus</i> (1)
<i>Campylobacter jejuni</i> (1)	<i>Salmonella dublin</i> (1)
<i>Candida albicans</i> (1)	<i>Salmonella typhimurium</i> (1)
<i>Citrobacter freundii</i> (1)	<i>Salmonella urbana</i> (1)
<i>Clostridium difficile</i> (2)	<i>Serratia liquefaciens</i> (3)
<i>Clostridium perfringens</i> (1)	<i>Shigella dysenteriae</i> (1)
<i>Enterobacter cloacae</i> (1)	<i>Shigella flexneri</i> (1)
<i>Enterobacter agglomerans</i> (1)	<i>Shigella sonnei</i> (1)
<i>Escherichia coli</i> non-O157:H7 (13)	<i>Staphylococcus aureus</i> (1)
<i>Escherichia fergusonii</i> (1)	<i>Staphylococcus aureus</i> (Cowan I) (1)
<i>Escherichia hermannii</i> (1)	<i>Staphylococcus epidermidis</i> (1)
<i>Hafnia alvei</i> (1)	<i>Streptococcus faecalis</i> (1)
<i>Helicobacter pylori</i> (1)	<i>Xanthomonas maltophilia</i> (1)
<i>Klebsiella pneumoniae</i> (1)	<i>Yersinia enterocolitica</i> (1)
<i>Providencia stuartii</i> (1)	
<i>Proteus vulgaris</i> (1)	

Three crossreactive bacteria have been identified: *Salmonella hilversum*, an *E. coli* O157:H38 and an *E. coli* O157:H1. *Salmonella hilversum* is not found in humans. *E. coli* O157:H38 and H1 are rare human isolates.

Twenty-one known laboratory *E. coli* O157:H7 isolates were shown to be correctly positive with the Immunocard STAT! *E. coli* O157 Plus test. No *E. coli* O157:H7 isolates have been observed to be negative.

Clinical testing with stool showed that Immunocard STAT! *E. coli* O157 Plus did not cross-react with a variety of other stool pathogens. A list of pathogens contained in stools which did not cross-react with the Immunocard test is given below:

Non Crossreacting Stool Pathogens (number of Stools)	
<i>Aeromonas</i> (1)	<i>E. coli</i> O157:H19 (toxin negative) (1)
<i>Campylobacter</i> (34)	<i>E. coli</i> O18:H7 (1)
<i>E. coli</i> O103:H2 (1)	<i>E. coli</i> O26:H11 (4)
<i>E. coli</i> O103:NM (1)	<i>E. coli</i> O69:H47 (1)
<i>E. coli</i> O111:NM (1)	<i>Plesiomonas shigelloides</i> (1)
<i>E. coli</i> O113:H21 (1)	<i>Salmonella</i> spp. (18)
<i>E. coli</i> O121:NM (1)	<i>Shigella</i> (7)
<i>E. coli</i> O145:NM (1)	<i>Yersinia</i> (1)

The following strains of non-Shiga toxin-producing *E. coli* O157 and non-O157 *E. coli* were spiked into positive and negative stool and tested with the Immunocard STAT! *E. coli* O157 Plus device at ≥ 6 x 10⁸ cells/mL. All were found to neither crossreact nor interfere with a positive result.

<i>E. coli</i> O157:H44	<i>E. coli</i> O26:H11
<i>E. coli</i> O157:H3	<i>E. coli</i> O111:NM
<i>E. coli</i> O157:H12	<i>E. coli</i> O125:NM
<i>E. coli</i> O126:H27	<i>E. coli</i> O55:B5:H-
<i>E. coli</i> O45:H2	<i>E. coli</i> #8739
<i>E. coli</i> ETEC LT	<i>E. coli</i> #9637
<i>E. coli</i> O121:H19	

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REFERENCES

1. Griffin PM, Ostroff SM, Tauxe RV, Green KD, Wells JG, Lewis JH et al. Illnesses associated with *Escherichia coli* O157:H7 infections. A broad clinical spectrum. Ann Intern Med. 1998; 109:705-712.
2. Karmali MB. Laboratory diagnosis of verotoxin-producing *Escherichia coli* infection. Clin Micro Newsletter. 1987; 9:65-70.
3. O'Brien AD, Holmes RK. Shiga and shiga-like toxins. Micro Reviews. 1987; 51:206-219.
4. Tarr P. *Escherichia coli* O157:H7: Clinical, diagnostic and epidemiological aspects of human infection. Clin Inf Dis. 1995; 20:1-10.
5. Alison O'Brien and James Kaper, eds., O'Brien AD, Kaper J. Shiga toxin-producing *Escherichia coli*: Yesterday, today and tomorrow. ASM, Washington, DC. 1998.
6. Nataro J, Kaper J. Diarrheagenic *Escherichia coli*. Clin Micro Rev. 1998; 11:142-201.
7. Itoh Y. et. al. Enterohemorrhagic *Escherichia coli* O157:H7 present in radish sprouts. Applied and Env. Microbiol. 1998; 64:1532-1535.
8. Chers ML, et al. An outbreak of *Escherichia coli* O157:H7 infection associated with leaf lettuce consumption. J Inf Dis. 1998; 177:1588-1593.
9. Manley WA, Walford JR, Harrison JF, McKee GL. ImmunoCard STAT! (ICS) *E. coli* O157:H7 test in a suspect waterborne outbreak in Wyoming. ASM, Abstract submitted for 1999 General Meeting. 1998.
10. Lembock AR, Willis DH. Prospective evaluation of a 10 minute ImmunoCard STAT! (ICS) *E. coli* O157:H7 test at a children's hospital. AMS, Abstract submitted for 1999 General Meeting. 1998.
11. Kay BA, Griffin PM, Stockbine VA, Wells JG. Too fast food: Bloody diarrhea and death from *Escherichia coli* O157:H7. Clin Micro Newsletter. 1994; 16:17-19.
12. Neill MA. *Escherichia coli* O157:H7. A pathogen of no small renown. Infec Dis Newsletter. 1991; 10:19-24.
13. Barrett B, Laurie KL, Pricou NJ. Evaluation of ImmunoCard STAT! STEC O157 for detection of *E. coli* O157:H7 in human stool specimens, enrichment broth and non-sorbitol fermenting *E. coli* from SMAC agar. ASM, Abstract submitted for 1999 General Meeting. 1998.
14. Stapp JR, Jelacic S, Fischer M, Yea YL, Clausen CR, Tarr PI. Sensitivity assessment of the Meridian ImmunoCard STAT! for the detection of *Escherichia coli* O157:H7 in children's stools. ASM, Abstract submitted for 1999 General Meeting. 1998.
15. Park CH, Alverson J. Evaluation of the Meridian ImmunoCard STAT! STEC O157 for the rapid detection of antigen directly from human stools and enhanced broth. ASM, Abstract submitted for 1999 General Meeting. 1998.
16. Mackenzie A, Chan F, Orrbine E, Hoban D, Kennedy W, Willis D, et al. and the CPKDR- co-investigators. Sensitivity and specificity of ImmunoCard STAT test for the detection of enterohemorrhagic *Escherichia coli* O157:H7 in stools. ASM, Abstract submitted for 1999 General Meeting. 1998.
17. Data on file, Meridian Bioscience, Inc.

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.