

Human CA15-3 LCR Calib. Grade

Cat. No. A32231H

STORE AT -30 to -10°C

DESCRIPTION Human Breast Tumor Antigen CA15-3 Calibrator Grade, Low Cross-Reactivity

TYPE Antigen

SHIPPING & STORAGE INFORMATION Store at -30 to -10°C. Avoid multiple freeze/thaw cycles.

SAFETY PRECAUTIONS Refer to the appropriate Safety Data Sheet (SDS) for additional information.

NOTES This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE Cell Culture Supernatant
FORMAT Calibrator Grade, Liquid
BUFFER Phosphate Buffered Saline, pH 7.4 ± 0.1
PRESERVATIVES 0.1% Sodium Azide
PURIFICATION Product is 0.2 µm filtered
INACTIVATION Not Applicable

Applications

Suitable for use in Chemiluminescent Immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.