

Human CA19-9 Antigen Grade

Cat. No. A86199H

STORE AT -30 to -10°C

DESCRIPTION	Human Cancer Antigen 19-9 (CA19-9) (Gastrointestinal Cancer) Antigen Grade
TYPE	Antigen
SHIPPING & STORAGE INFORMATION	Store at -30 to -10°C. Avoid multiple freeze/thaw cycles.
SAFETY PRECAUTIONS	Refer to the appropriate Safety Data Sheet (SDS) for additional information.
NOTES	This reagent has been manufactured under ISO13485 quality management system and is suitable for research or further manufacturing use only. Centrifuge before opening to ensure complete recovery of vial contents.

» For ordering and to learn about pack sizes, please [click here](#).

General Features/Characteristics

SOURCE	Human metastatic liver carcinoma
FORMAT	Purified, Liquid
PURIFICATION	Gel filtration and column chromatography
BUFFER	Phosphate Buffered Saline, pH 7.4 containing 3% Sucrose
PRESERVATIVES	0.05% Sodium Azide
INACTIVATION	Not Applicable

Applications

Immunogen, Calibrator in Immunoassays and in vitro studies. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded